

Purified Human Pancreatic Islets, CIT Purification Density Gradients – *A Standard Operating Procedure of the NIH Clinical Islet Transplantation Consortium*

The NIH CIT Consortium Chemistry Manufacturing Controls Monitoring

Committee: J. Ansite, A.N. Balamurugan, B. Barbaro, J. Battle, D. Brandhorst, J. Cano, X. Chen, S. Deng, D. Feddersen, A. Friberg, T. Gilmore, J.S. Goldstein, E. Holbrook, A. Khan, T. Kin, J. Lei, E. Linetsky, C. Liu, X. Luo, K. McElvaney, Z. Min, J. Moreno, D. O’Gorman, K.K. Papas, G. Putz, C. Ricordi, G. Szot, T. Templeton, L. Wang, J.J. Wilhelm, J. Willits, T. Wilson, X. Zhang

The NIH CIT Consortium

Emory University: J. Avila, B. Begley, J. Cano, S. Carpentier, E. Holbrook, J. Hutchinson, C.P. Larsen, J. Moreno, M. Sears, N.A. Turgeon, D. Webster

Massachusetts General Hospital: S. Deng, J. Lei, J.F. Markmann

NIAID: N.D. Bridges, C.W. Czarniecki, J.S. Goldstein, G. Putz, T. Templeton, T. Wilson

NIDDK: T.L. Eggerman

Northwestern University: P. Al-saden, J. Battle, X. Chen, A. Hecyk, H. Kissler, X. Luo, M. Molitch, N. Monson, E. Stuart, A. Wallia, L. Wang, S. Wang, X. Zhang

University of Alberta, Edmonton: D. Bigam, P. Campbell, P. Dinyari, T. Kin, N. Kneteman, J. Lyon, A. Malcolm, D. O’Gorman, C. Onderka, R. Owen, R. Pawlick, B. Richer, S. Rosichuk, D. Sarman, A. Schroeder, P.A. Senior, A.M.J. Shapiro, L. Toth, V. Toth, W. Zhai

University of California–San Francisco: K. Johnson, J. McElroy, A.M. Posselt, M. Ramos, T. Rojas, P.G. Stock, G. Szot

University of Illinois, Chicago: B. Barbaro, J. Martellotto, J. Oberholzer, M. Qi, Y. Wang

University of Iowa (Data Coordinating Center): L. Bayman, K. Chaloner, W. Clarke, J.S. Dillon, C. Diltz, G.C. Doelle, D. Ecklund, D. Feddersen, E. Foster, L. G. Hunsicker, C. Jasperson, D-E Lafontant, K. McElvaney, T. Neill-Hudson, D. Nollen, J. Qidwai, H. Riss, T. Schwieger, J. Willits, J. Yankey

University of Miami: R. Alejandro, A.C. Corrales, R. Faradji, T. Froud, A.A. Garcia, E. Herrada, H. Ichii, L. Inverardi, N. Kenyon, A. Khan, E. Linetsky, J. Montelongo, E. Peixoto, K. Peterson, C. Ricordi, J. Szust, X. Wang

University of Minnesota: M.H. Abdulla, J. Ansite, A.N. Balamurugan, M.D. Bellin, M. Brandenburg, T. Gilmore, J. V. Harmon, B.J. Hering, R. Kandaswamy, G. Loganathan, K. Mueller, K.K. Papas, J. Pedersen, J.J. Wilhelm, J. Witson

University of Pennsylvania: C. Dalton-Bakes, H. Fu, M. Kamoun, J. Kearns, Y. Li, C. Liu, E. Luning Prak, Y. Luo, E. Markmann, Z. Min, A. Naji, M. Palanjan, M. Rickels, R. Shlansky-Goldberg, K. Vivek, A.S. Ziaie

University of Wisconsin: L. Fernandez, D.B. Kaufman, L. Zitur

Uppsala University: D. Brandhorst, A. Friberg, O. Korsgren

Supported by grants from the National Institute of Allergy and Infectious Diseases and the National Institute for Diabetes and Digestive and Kidney Diseases.

- At Emory University, U01AI089317.
- At Northwestern University, U01AI089316.
- At the University of Alberta, Edmonton: U01AI065191.
- At the University of California, San Francisco, U01DK085531.
- At the University of Illinois, Chicago, 5U01DK070431-10.
- At the University of Iowa, U01DK070431.
- At the University of Miami, U01DK070460.
- At the University of Minnesota, U01AI065193.
- At the University of Pennsylvania, U01DK070430.
- At Uppsala University, U01AI065192.

In addition, the study was supported by the following GCRC and CTSA awards:

- At Emory University: UL1TR000454.
- At Northwestern University: 5UL1RR025741 and 8UL1TR000150.
- At the University of California, San Francisco, UL1TR000004.
- At the University of Illinois, Chicago, UL1TR000050.
- At the University of Miami: 1UL1TR000460.
- At the University of Minnesota: 5M01-RR000400 and UL1TR000114.
- At the University of Pennsylvania: UL1TR000003.

Address correspondence to: Camillo Ricordi MD, Chairman, CIT Steering Committee, ricordi@miami.edu

To cite this article

Purified Human Pancreatic Islets, CIT Purification Density Gradients – A Standard Operating Procedure of the NIH Clinical Islet Transplantation Consortium

CellR4 2014; 2 (3): e991

DAIT, NIAID, NIH

**SOP
ATTACHMENT**



Document No. 3106, B10	Revision No. 00	Effective Date 07 July 2008	Supersedes Date N/A	Page 1 of 2
Document Title: PURIFIED HUMAN PANCREATIC ISLETS, CIT PURIFICATION DENSITY GRADIENTS				

Manufacturing Site: _____ Date: _____

1. Materials:

Material	Source	For #1, Preparation Date For #2 & #3, Lot #	Expiration Date
1. CIT Purification Solution	Made In-house		
2. OptiPrep	Nycomed, Product #1114542		
3. Gradient Stock Solution	Mediatech, Product No. 99-674-CM		

2. Procedure

2.1 **High Purification Density Gradient, 1.10 g/mL ± 0.01 g/mL**

Material	Quantity Required	Initial Quantity Used	Additional Quantity Used
CIT Purification Solution	81.25 g	g	g
Gradient Stock Solution	468.75 g	g	g

- 2.1.1 In a BSC place a sterile 500 mL bottle on a balance.
- 2.1.2 Add 81.25 g of CIT Purification Solution to the bottle.
- 2.1.3 Add 468.75 g of Gradient Stock Solution to the bottle.
- 2.1.4 Cap the bottle and mix by gentle inversion at least five times.
- 2.1.5 Check that the weight is 550.0 g. Observed Weight: _____ g
- 2.1.6 Check the density of the solution with a densitometer.
Measured Density: _____ g/mL
If the density is not 1.09 g/mL to 1.11 g/mL, the solution's density may be adjusted by adding the higher or lower density component. Record the weight of component added in the table, above. If the density cannot be adjusted into the range, the solution may not be used. Initiate and document an investigation of the process and materials.

Islets Lot Number: _____

Document No. 3106, B10	Revision Number 00	Effective Date 07 July 2008	Supersedes Date N/A	Page 2 of 2
Document Title: PURIFIED HUMAN PANCREATIC ISLETS, CIT PURIFICATION DENSITY GRADIENTS				

- 2.1.7 Label the bottle with:
- “High Purification Density Gradient, 1.10 g/mL”
 - Islets Lot Number
 - “Store at 2°C to 8°C”
 - Date Prepared
 - Expiration Date (the end of the day after preparation)
 - Initials of the person who prepared the gradient
- 2.1.8 Store the bottle of solution at 2°C to 8°C before use.

Total volume prepared: _____ mL

Prepared by: _____ Date: _____

Reviewed by: _____ Date: _____

2.2 Low Purification Density Gradient **1.06 g/mL** ± 0.01 g/mL

Material	Quantity Required	Initial Quantity Used	Additional Quantity Used
CIT Purification Solution	245.7 g	g	g
OptiPrep	29.9 g	g	g

- 2.2.1 In a BSC place a sterile 500 mL bottle on a balance.
- 2.2.2 Add 245.7 g of CIT Purification Solution to the bottle.
- 2.2.3 Add 29.9 g of OptiPrep to the bottle.
- 2.2.4 Cap the bottle and mix by gentle inversion at least five times.
- 2.2.5 Check that the weight is 275.6 g. Observed Weight: _____ g
- 2.2.6 Check the density of the solution with a densitometer.
Measured Density: _____ g/mL
If the density is not 1.05 g/mL to 1.07 g/mL, the solution’s density may be adjusted by adding the higher or lower density component. Record the weight of component added in the table, above. If the density cannot be adjusted into the range, the solution may not be used. Initiate and document an investigation of the process and materials.
- 2.2.7 Label the bottle with:
- “Low Purification Density Gradient, 1.06 g/mL”
 - Islets Lot Number
 - “Store at 2°C to 8°C”
 - Date Prepared
 - Expiration Date (the end of the day after preparation)
 - Initials of the person who prepared the gradient
- 2.2.8 Store the bottle of solution at 2°C to 8°C before use.

Total volume prepared: _____ mL

Prepared by: _____ Date: _____

Reviewed by: _____ Date: _____

Islets Lot Number: _____