

The NIH CIT Consortium Chemistry Manufacturing Controls Monitoring Committee:

J. Ansite, A.N. Balamurugan, B. Barbaro, J. Battle, D. Brandhorst, J. Cano, X. Chen, S. Deng, D. Feddersen, A. Friberg, T. Gilmore, J.S. Goldstein, E. Holbrook, A. Khan, T. Kin, J. Lei, E. Linetsky, C. Liu, X. Luo, K. McElvaney, Z. Min, J. Moreno, D. O’Gorman, K.K. Papas, G. Putz, C. Ricordi, G. Szot, T. Templeton, L. Wang, J.J. Wilhelm, J. Willits, T. Wilson, X. Zhang

The NIH CIT Consortium

Emory University: J. Avila, B. Begley, J. Cano, S. Carpentier, E. Holbrook, J. Hutchinson, C.P. Larsen, J. Moreno, M. Sears, N.A. Turgeon, D. Webster

Massachusetts General Hospital: S. Deng, J. Lei, J.F. Markmann

NIAID: N.D. Bridges, C.W. Czarniecki, J.S. Goldstein, G. Putz, T. Templeton, T. Wilson

NIDDK: T.L. Eggerman

Northwestern University: P. Al-saden, J. Battle, X. Chen, A. Hecyk, H. Kissler, X. Luo, M. Molitch, N. Monson, E. Stuart, A. Wallia, L. Wang, S. Wang, X. Zhang

University of Alberta, Edmonton: D. Bigam, P. Campbell, P. Dinyari, T. Kin, N. Kneteman, J. Lyon, A. Malcolm, D. O’Gorman, C. Onderka, R. Owen, R. Pawlick, B. Richer, S. Rosichuk, D. Sarman, A. Schroeder, P.A. Senior, A.M.J. Shapiro, L. Toth, V. Toth, W. Zhai

University of California–San Francisco: K. Johnson, J. McElroy, A.M. Posselt, M. Ramos, T. Rojas, P.G. Stock, G. Szot

University of Illinois, Chicago: B. Barbaro, J. Martellotto, J. Oberholzer, M. Qi, Y. Wang

University of Iowa (Data Coordinating Center): L. Bayman, K. Chaloner, W. Clarke, J.S. Dillon, C. Diltz, G.C. Doelle, D. Ecklund, D. Feddersen, E. Foster, L. G. Hunsicker, C. Jaspersen, D-E Lafontant, K. McElvaney, T. Neill-Hudson, D. Nollen, J. Qidwai, H. Riss, T. Schwieger, J. Willits, J. Yankey

University of Miami: R. Alejandro, A.C. Corrales, R. Faradji, T. Froud, A.A. Garcia, E. Herrada, H. Ichii, L. Inverardi, N. Kenyon, A. Khan, E. Linetsky, J. Montelongo, E. Peixoto, K. Peterson, C. Ricordi, J. Szust, X. Wang

University of Minnesota: M.H. Abdulla, J. Ansite, A.N. Balamurugan, M.D. Bellin, M. Brandenburg, T. Gilmore, J. V. Harmon, B.J. Hering, R. Kandaswamy, G. Loganathan, K. Mueller, K.K. Papas, J. Pedersen, J.J. Wilhelm, J. Witson

University of Pennsylvania: C. Dalton-Bakes, H. Fu, M. Kamoun, J. Kearns, Y. Li, C. Liu, E. Luning-Prak, Y. Luo, E. Markmann, Z. Min, A. Naji, M. Palanjan, M. Rickels, R. Shlansky-Goldberg, K. Vivek, A.S. Ziaie

University of Wisconsin: L. Fernandez, D.B. Kaufman, L. Zitur

Uppsala University: D. Brandhorst, A. Friberg, O. Korsgren

Supported by grants from the National Institute of Allergy and Infectious Diseases and the National Institute for Diabetes and Digestive and Kidney Diseases.

- At Emory University, U01AI089317.
- At Northwestern University, U01AI089316.
- At the University of Alberta, Edmonton: U01AI065191.
- At the University of California, San Francisco, U01DK085531.
- At the University of Illinois, Chicago, 5U01DK070431-10.
- At the University of Iowa, U01DK070431.
- At the University of Miami, U01DK070460.
- At the University of Minnesota, U01AI065193.
- At the University of Pennsylvania, U01DK070430.
- At Uppsala University, U01AI065192.

In addition, the study was supported by the following GCRC and CTSA awards:

- At Emory University: UL1TR000454.
- At Northwestern University: 5UL1RR025741 and 8UL1TR000150.
- At the University of California, San Francisco, UL1TR000004.
- At the University of Illinois, Chicago, UL1TR000050.
- At the University of Miami: 1UL1TR000460.
- At the University of Minnesota: 5M01-RR000400 and UL1TR000114.
- At the University of Pennsylvania: UL1TR000003.

Address correspondence to: Camillo Ricordi MD, Chairman, CIT Steering Committee,
ricordi@miami.edu

To cite this article

Purified Human Pancreatic Islets Master Production Batch Record, Part 1 University of Illinois, Chicago & University of Miami (Product Codes PHPI-A-01 & PHPI-L-01) – Standard Operating Procedure of the NIH Clinical Islet Transplantation Consortium

CellR4 2017; 5 (2): e2286

Document No.
SOP 3101, B02-1Revision No.
07Effective Date
06 Aug 2011Supersedes Date
02 May 2011

Page 1 of 39

Document Title:

**PURIFIED HUMAN PANCREATIC ISLETS
MASTER PRODUCTION BATCH RECORD, PART 1
UNIVERSITY OF ILLINOIS, CHICAGO & UNIVERSITY OF MIAMI
(PRODUCT CODES PHPI-A-01 & PHPI-L-01)**

1.0 MASTER PRODUCTION BATCH RECORD APPROVAL*Approvals on file*

Camillo Ricordi, M.D.
University of Miami, Miami, Florida

Date: _____

James F. Markmann, M.D., Ph.D.
Massachusetts General Hospital, Boston, Massachusetts

Date: _____

Jose Oberholzer, M.D.
University of Illinois at Chicago

Date: _____

Christine W. Czamiecki, Ph.D.
DAIT, NIAID, NIH, Bethesda, Maryland

Date: _____


Changes to this Master Production Batch Record must be proposed to the Chief, Regulatory Affairs, DAIT, NIAID, NIH, and approved by all the original signatories, or their successors, before implementation.

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 2 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				


2.0 FLOWCHART AND SAMPLING TABLE

2.1 Production Process Flowchart (MPBR)




Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 3 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				



Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 4 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				



Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 5 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

2.2 Samples and Tests

MPBR SECTION	SAMPLE TYPES & QUANTITIES	
	PROCESS CONTROL TESTS	TESTS
5.1	Preservation Solution, ≥ 3 mL	Sterility (21 CFR 610.12) & Fungal Culture
7.1.3	Pancreas Digest, ≤ 1 -2 mL periodically	Acinar Amount, # of Islets, % Free Islets, % Fragmented
7.5.1	Diluted Pancreas Digest, 2 X 100 μ L	Islets Count
8.3.7	Purification Fractions, 0.5 mL/each of 12 fractions & 0.5 mL of W1 fraction, each COBE Run	Islets Purity (%)
8.4.3	Supplementary Purification Islets, 2 X 100 μ L (Optional)	Islets Count
10.2	Purified Islets, 2 X 100 μ L, High, Middle, Low Purity Levels	Islets Count
12.10	Cultured Islets, All Measured, High, Middle, Low Purity Levels	Settled Tissue Volume
12.13	Cultured Islets, 2 X 100 μ L, High, Middle, Low Purity Levels	Post-culture Islets Count
	INTERIM & FINAL CERTIFICATES OF ANALYSIS	
11.1	Suspension, 400 IEQ, High Purity Islets	Glucose Stimulated Insulin Release
12.11.5	Supernatant above cultured islets, volume according to institution's procedure, High, Middle, Low Purity Levels	Gram Stain
12.13 & 12.14, or 12.17.1	Suspension, 2 X 100 μ L/Each Final Product T-75 Flask	Islets Identity, Quantity, Concentration
12.17.2	Suspension, 100 IEQ/Each Final Product T-75 Flask	Viability
12.17.3	Supernatant above cultured islets, 1 mL/Each Final Product T-75 Flask	Endotoxin
12.18	Combined Islets, All Measured, High, Middle, Low Purity Levels	Settled Tissue Volume
	FINAL CERTIFICATE OF ANALYSIS ONLY	
12.14	Suspension, 400 IEQ, High Purity Islets (Post-culture Sample)	Glucose Stimulated Insulin Release
12.17.2	Volume according to institution's procedure of islets suspension in each T-75 Flask	Sterility (21 CFR 610.12) & Fungal Culture
	REQUIRED PRODUCT CHARACTERIZATION TESTS FOR INFORMATION ONLY	
5.6	Superficial biopsy of approximately 3 mm X 3 mm X 3 mm	MCP-1 and Tissue Factor
12.14	Suspension, 4,000 IEQ, High Purity Islets	<i>In vivo</i> (Nude Mouse) Islets Function
12.17.2	Suspension, 1,000 IEQ/Each Final Product T-75 Flask	Cell Composition
12.17.2	Suspension, 500 to 1,000 IEQ/Each Final Product T-75 Flask	MCP-1 and Tissue Factor
12.17.2	Suspension, 4 X 500 IEQ from T-75 Flask #1 in 1.8 mL cryovials	NIDDK Repository
	OPTIONAL PRODUCT CHARACTERIZATION TESTS FOR INFORMATION ONLY	
11.1	Suspension, 3 X 100 IEQ, High Purity Islets	Pre-culture DNA Content
11.1	Suspension, 3 X 100 IEQ, High Purity Islets	Nuclei Measurement
12.14	Suspension, 3 X 100 IEQ, High Purity Islets	Post-culture DNA Content
12.14	Suspension, 3 X 100 IEQ, High Purity Islets	Nuclei Measurement
12.14	Suspension, 500 IEQ, High Purity Islets	ATP/DNA
12.14	Suspension, 5,000 IEQ, High Purity Islets	OCR/DNA
12.14	Suspension, 5,000 IEQ, High Purity Islets	Molecular Profiling
12.14	Suspension, 500 IEQ, High Purity Islets	Islets Fraction
12.17.2	Suspension, 2,000 IEQ/Each Final Product T-75 Flask	β -cell Viability

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 6 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Note: Materials used in this process may transmit infectious agents. Therefore, each person participating in this process must be trained in, and follow, the institution's procedures for handling potentially infectious agents. All waste materials from this process that may have contacted the pancreas or the islets must be discarded as Biohazardous Waste.

Note: It is extremely important to protect the pancreas and the islets from contamination by adventitious microorganisms and pyrogenic agents. Reagents and equipment that may contact the pancreas or islets must be sterile, pyrogen-free, and single-use whenever possible. The institution's procedures for aseptic technique must be followed throughout the execution of this Production Batch Record. All "open" procedure steps must be performed in a clean and disinfected Certified Class II area or Biological Safety Cabinet (BSC).

Note If, at any time during the execution of this Production Batch Record, you observe:

- 1) potential discrepancies in the identification of the pancreas or islets,*
- 2) unusual appearance of any materials,*
- 3) unusual, or improper performance of any equipment, or*
- 4) inadvertent deviations from the process as defined in this Production Batch Record or the institution's established procedures;*

you must notify the Laboratory Director, or designee, immediately.

The Laboratory Director, or designee, must investigate the observation, and write, sign and date a report giving the details of the observation and its resolution according to the institution's procedures. The occurrence of the event is documented in this Production Batch Record by writing "See Report #X" at the location in the Batch Record where the observation occurred. When allowed by the institution's procedures the report, or a copy, must be filed with this Batch Record. When not allowed, it must be traceable through the unique identification number ("Report #X") written in the Batch Record. The process for reporting a deviation to the CCMC as defined in DAIT SOP 3200 must also be followed.

3.0 LABORATORY PREPARATION

- 3.1 Identification of Institution, Personnel, Raw Materials and Purchased Reagents, Sterilized Items, Equipment and Disposable Items

- 3.1.1 Institution Manufacturing Purified Human Pancreatic Islets Product

Name of Institution: _____

- 3.1.2 Personnel

Attach to this Batch Record a list of the names of all personnel directly involved in the execution of this Batch Record and their signatures and initials, or have them sign and initial the table below.

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 7 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

PRINTED NAME	SIGNATURE	INITIALS

3.1.3 Raw Materials and Purchased Reagents

Below is a list of the raw materials and purchased reagents used in this procedure, including their catalog numbers and suppliers, where specific Catalog Numbers and Suppliers are required. Record in the table the Catalog Number and Supplier, where not already specified, and the lot number and expiration date of each material used.

RAW MATERIAL AND PURCHASED REAGENTS	CATALOG NUMBER	SUPPLIER	LOT NUMBER	EXPIRATION DATE
1. CMRL 1066, Supplemented, CIT Modifications				
2. CMRL 1066 Transplant Media, contains Hepes and without Sodium Bicarbonate				
3. Hanks' Balanced Salt Solution (HBSS), 1X				
4. Heparin Sodium Injection USP, Preservative Free		_____ Units/mL		
5. HEPES Buffer, 1 M				
6. Gradient Stock Solution				
7. Phase I Solution				
8. Cold Storage/Purification Stock Solution				
9. Albumin Human USP, 25% Solution				
10. Hydrochloric Acid NF, 1 N				
11. Insulin-like Growth Factor-1 (IGF-1), 1.0 mg/vial	CM001	Cell Sciences		
12. Insulin Human Injection USP, Recombinant				

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 8 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

RAW MATERIALS AND PURCHASED REAGENTS (Continued)

RAW MATERIAL AND PURCHASED REAGENTS	CATALOG NUMBER	SUPPLIER	LOT NUMBER	EXPIRATION DATE
13a. Collagenase NB 1 GMP Grade	N0002937	SERVA/Nordmark		
13b. Neutral Protease NB GMP Grade	N0002936	SERVA/Nordmark		
14a. Collagenase NB 1 Premium Grade	17455	SERVA/Nordmark		
14b. Neutral Protease NB	30301	SERVA/Nordmark		
15a. Clzyme Collagenase HA	001-1000	VitaCyte LLC		
15b. Clzyme Thermolysin	002-1000	VitaCyte LLC		
16. Liberase MTF C/T GMP Grade	05339880001	Roche Diagnostics		
17. OptiPrep				
18. Trimming Solution				
19. Human Pancreas, Deceased Donor	See Section 4.2 and SOP 3108			
20. PentaStarch, 10% Solution				
21. Povidone Iodine USP, 10%				
22. Pulmozyme (dornase alpha), 2.5 mL/vial, 1 mg/mL	NDC No. 50242-100-40	Genentech		
23. RPMI 1640 with L-Glutamine				
24. Sterile Water for Injection USP				
25. Viaspan (UW Solution)				
26. Biocoll Separating Solution, Density 1.100	L6155	Biochrome AG/ Cedarlane		
27. Stock Polysucrose Solution, sterile	99-662-CVS	Mediatech		
28. Islet Gradient 1.037, sterile	99-690-CIS	Mediatech		
29. Islet Gradient 1.096, sterile	99-691-CIS	Mediatech		
30. Islet Gradient 1.108, sterile	99-692-CIS	Mediatech		
31. Calcium Chloride USP (Dihydrate) (CaCl ₂ 2 H ₂ O)				
32. Calcium Chloride Injection USP				
33. Cefazolin Sodium USP				
34. Lisofylline, 60 mg/mL	Formula #0109-00	DiaKine Therapeutics		
35. Infusion Bag				

Verified by: _____

Date: _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 9 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

3.1.4 Sterilized Items

Attach a list of all items used in this process that have been sterilized, the sterilizer load numbers and dates, and verify that the sterilizations were performed within the time period validated by the institution.

Verified by: _____ **Date:** _____

3.1.5 Equipment

Attach a list of all equipment used in the manufacturing process, including identification numbers, serial numbers, etc.

Verified by: _____ **Date:** _____

3.1.6 Disposable Items

Attach a list of all disposable items used in this process, the supplier of each, the lot number, and the expiration date.

Verified by: _____ **Date:** _____

3.2 Biological Safety Cabinet and Laboratory Preparation

Prepare the laboratory, including the Biological Safety Cabinet (BSC), for islet isolation according to the institution's procedure(s) and record the preparation on the appropriate form(s) or logbook(s). Submit copies of the form(s) or logbook page(s) with this Batch Record.

Verified by: _____ **Date:** _____

3.3 Dilution Media Preparation

3.3.1 Equilibrate RPMI 1640 for digest dilution to room temperature prior to use for approximately 1 to 2 hours.

3.3.2 Prepare four 1L containers ahead of time and store at 2°C to 8°C before use:

REQUIRED	USED
1st Container	
400 mL of RPMI 1640	mL
200 mL of Albumin Human USP, 25% Solution	mL
200 Units of insulin (final concentration: 0.2 Units/mL)	Units
10,000 Units of heparin (final concentration: 10 Units/mL)	Units
2nd Container	
400 mL of RPMI 1640	mL
200 mL of Albumin Human USP, 25% Solution	mL
200 Units of insulin (final concentration: 0.2 Units/mL)	Units
10,000 Units of heparin (final concentration: 10 Units/mL)	Units

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 10 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

3rd Container	
500 mL of RPMI 1640	mL
100 mL of Albumin Human USP, 25% Solution	mL
200 Units of insulin (final concentration: 0.2 Units/mL)	Units
10,000 Units of heparin (final concentration: 10 Units/mL)	Units
4th Container	
500 mL of RPMI 1640	mL
100 mL of Albumin Human USP, 25% Solution	mL
200 Units of insulin (final concentration: 0.2 Units/mL)	Units
10,000 Units of heparin (final concentration: 10 Units/mL)	Units

Performed by: _____ **Date:** _____

Verified by: _____ **Date:** _____

- 3.3.3 Fill as many additional containers as needed with enough Albumin Human USP, 25% Solution each to provide a final concentration of 1.5% Albumin.

Number of additional containers: _____

Volume of each additional container: _____ mL

Volume collected in each additional container: _____ mL

Volume of Albumin Human USP, 25% Solution in each additional container _____ mL

Performed by: _____ **Date:** _____

Verified by: _____ **Date:** _____

4.0 PANCREAS ACCEPTANCE AND RECEIPT

- 4.1 Time of pancreas receipt in the lab: _____ (Record all times using the 24-hour clock)

Received by: _____ **Date:** _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 11 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

4.2 Pancreas Donor Qualification Record (NA = Not Available)

REQUIREMENTS A qualified donor must have "Yes" responses to all of the Inclusion Criteria (A), and "No" responses to all of the Exclusion Criteria (B & C).			
	Yes	No	NA
Container Label must specify Human Pancreas, and a UNOS or DDD number must be present.			
The Organ Procurement Organization (OPO) must be identified.			
A. Inclusion Criteria (The donor or pancreas must meet these criteria.)			
1. Pancreas Preservation in (i) UW, (ii) PF/UW, (iii) HTK, or (iv) PF/HTK Solution(s)			
2. Maximum 12 hour cold ischemia time			
3. Donor age 15-65 years			
4. Cause and circumstances of death acceptable to the transplant team			
B. Exclusion Criteria (Is there evidence of the following conditions?)			
1. History or biochemical evidence of Diabetes mellitus Type 1 or 2 (Transplant teams may consider donor HbA1C > 6.1% in the absence of transfusions in the week prior to death as an indication for exclusion, with discretion for donors who have received transfusions.)			
2. Pancreas from non-heart-beating cardiac death donors.			
3. Malignancies, other than resected basal squamous cell carcinoma or intracranial tumor as the cause of death			
4. Suspected or confirmed sepsis			
5. Evidence of clinical or active viral Hepatitis [A, B (HBcAg), C]. HBsAb+ is acceptable, if there is a history of vaccination.			
6. Acquired Immunodeficiency Syndrome (AIDS)			
7. HIV seropositivity (HIV-I or HIV-II), or HIV status unknown*			
8. HTLV-I or HTLV-II (Optional)			
9. Syphilis (RPR or VDRL positive)*			
10. Active viral encephalitis or encephalitis of unknown origin			
11. TSE or Creutzfeldt-Jacob Disease			
12. Suspected Rabies Diagnosis			
13. Treated or Active Tuberculosis			
14. Individuals who have received pit-hGH (pituitary growth hormone)			
15. Any medical condition that, in the opinion of the transplant team, precludes a reasonable possibility of a favorable outcome of the islet transplant procedure			
16. Clinical history and/or laboratory testing suggestive of West Nile Virus, Vaccinia, or SARS			
C. Exclusion Criteria – Behavioral Profiles (Is there evidence of the following conditions?)			
17. High-risk sexual behavior within 5 years prior to time of death: men who have had sex with men, individuals who have engaged in prostitution, and individuals whose sexual partners have engaged in high-risk sexual behavior			
18. Non-medical intravenous, intramuscular, or subcutaneous drug use within the past five years			
19. Persons with hemophilia or related clotting disorders who have received human-derived clotting factor concentrates			
20. Findings on history or physical examination consistent with an increased risk of HIV exposure			
21. Current inmates of correctional systems and individuals who have been incarcerated for more than 72 consecutive hours during the previous 12 months			

*Test results for Exclusion Criteria B. 7 and 9 are required by FDA regulation.

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 12 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Is donor qualified as pancreas source? Yes No (Circle One)

Recorded by: _____ **Date:** _____

Review by: _____ **Date:** _____

- 4.3 Examine the container in which the pancreas arrived and its label. Is the container clean, intact and labeled with the UNOS or DDD number that has been accepted and are a proper name and donor records present?

Yes No (Circle One)

Is the product packaged properly?

Yes No (Circle One)

Comments: _____

Examined by: _____ **Date:** _____

- 4.4 Record the following information from donor records provided by the OPO:

PANCREAS DONOR INFORMATION (NA = Not Available)

	OBSERVED	ACCEPTABLE?		
		Yes	No	NA
UNOS or DDD Number				
Name and Location of OPO				
OPO Unique Identifier (if applicable)				
Donor Consent for Islets Transplant Present				
Donor's Date of Birth				
Donor's Gender				
Donor's ABO				
Donor's Weight				
Donor's Height				
Donor's Body Mass Index				
Extent of Hemodilution (See Flowchart & Worksheet at the end of this document)				
Donor's CMV Status				

Recorded by: _____ **Date:** _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 13 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

5.0 PANCREAS PREPARATION

5.1 In-process Samples for Sterility Testing of Preservation Solution

Preservation Method: _____

Using sterile technique, open the pancreas container in a Class 100 area. Aseptically take at least a 3 mL sample of the preservation solution in which the pancreas was transported. Prepare and label the sample according to the institution's procedure and submit for sterility (21 CFR 610.12) and fungal culture testing to the appropriate laboratory. Attach a copy of the requisition form to the Production Batch Record.

Sample Collected by: _____ **Date:** _____

Record the test results, when available, in Section 17.1.

Note: In some cases pancreas cleaning and cannulation are partially or completely performed immediately after the pancreas is procured and before it is delivered to the lab. In these cases, records of these activities will be made and filed with this Production Batch Record.

5.2 Move the pancreas to a cold tray containing Trimming Solution plus 1 g/L Cefazolin Sodium USP and remove excess tissue.

Process Start time: _____

Performed by: _____ **Date:** _____

5.3 Examine the cleaned pancreas and record observations in the table below.

Check only one line in each category.

Fat	___ Clean	Edema	___ None
	___ Average		___ Interstitial Edema
	___ Patchy Infiltration		___ Slight Overall Swelling
	___ Heavily Infiltrated		___ Overly Distended
Flush	___ Well Flushed	Texture	___ Very Soft
	___ Poorly Flushed		___ Soft
			___ Firm (normal)
			___ Many Firm Areas (Fibrotic)
			___ Rigid Throughout
Blood	___ Blood on Capillaries	Pancreas Condition	___ Intact
	___ Blood in Intra-Parenchymal		___ Capsular Damage
	___ No Blood Present		___ Parenchymal Damage

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 14 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Gross pathology observed? Yes No (Circle One)

Comments: _____

Examined by: _____ **Date:** _____

- 5.4 Prepare the CIT Digestion Solution according to DAIT SOP 3106, B01, and file the record of preparation with this Batch Record.

Performed by: _____ **Date:** _____

- 5.5 Optional Pancreas Surface Decontamination

If desired, place the pancreas in 250 mL of HBSS or preservation solution containing 1 mg/mL Cefazolin Sodium USP, or in 250 mL of 10% Povidone Iodine USP solution. Rinse the pancreas with 400 mL of plain HBSS 1X, transfer it to a new container of 400 mL of plain HBSS 1X, and rinse again. Remove the original pan and instruments from the BSC, and replace with clean, sterile pan and instruments.

Pancreas surface decontamination method: _____

Documented by: _____ **Date:** _____

- 5.6 Pancreas Biopsy

Collect a superficial biopsy of approximately 3 mm X 3 mm X 3 mm from the area within 1 cm of the main duct of the donor pancreas for required product characterization MCP-1 and tissue factor testing. Prepare and ship the sample according to instructions in the CIT Islets Lab Binder. Report the results in PBR Section 17.3.

Performed by: _____ **Date:** _____

- 5.7 Pancreas Weight

After excess tissue is trimmed from the pancreas, weigh the pancreas.

Initial Trimmed Pancreas Weight: _____ g

Recorded by: _____ **Date:** _____

Verified by: _____ **Date:** _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 15 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

5.8 CIT Enzyme Solution Preparation

Prepare the CIT Enzyme Solution described in the appropriate procedure reference below. Cross out the two references not used.

5.8.1 Prepare the CIT Enzyme Solution – SERVA Enzymes according to DAIT SOP 3106, B11.

5.8.2 Prepare the CIT Enzyme Solution – Vitacyte Enzymes and VitaCyte/SERVA Enzymes Combination according to DAIT SOP 3106, B13.

5.8.3 Prepare the CIT Enzyme Solution – Roche Enzymes according to DAIT SOP 3106, B14.

File the record of CIT Enzyme Solution preparation with this Batch Record.

Recorded by: _____ **Date:** _____

5.9 CIT Enzyme Solution (Specify Units of each enzyme)

Collagenase Activity actually used: _____

Neutral Protease Activity actually used: _____

Thermolysin Activity actually used: _____

Cross out the line above not used.

CIT Enzyme Solution volume prepared: _____ mL

Verified by: _____ **Date:** _____

5.10 Pancreas Cannulation

The pancreas will be perfused in a controlled manner, using separate cannulae for the head and tail. After the pancreas is cleaned of excess tissue, cut the pancreas to separate the head and tail, and cannulate the main pancreatic ducts with 16 to 22 gauge cannulae, one at the head and one at the tail. You may use a small cannula as a thread down the duct from the head of the pancreas to facilitate the identification of the duct for the cannulation process.

Performed by: _____ **Date:** _____

5.11 After the two portions of pancreas are cannulated, continue to remove excess tissue if necessary. Place this additional trimmed tissue in a tared container.

Comments on pancreas receipt and preparation for perfusion: _____

Written by: _____ **Date:** _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 16 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

6.0 PANCREAS PERFUSION

6.1 Assemble perfusion equipment according to the institution's procedure.

Performed by: _____ **Date:** _____

6.2 Perfuse the pancreas with the CIT Enzyme Solution.

- If indicated by the institution's procedures, prime the perfusion circuit by pumping HBSS, 1X, through it. Confirm the absence of leaks or loose connections, and drain the perfusion circuit.
- Add CIT Enzyme Solution (Section 5.5) at 4°C to 8°C to the chamber and refill the perfusion circuit with it. Remove all air bubbles.
- Connect the perfusion tubing to the cannula and perfuse the pancreas for 4 to 10 minutes at 60 to 80 mm Hg, followed by 4 to 6 minutes (8 minutes maximum in case of poor distension) at 160 to 180 mm Hg at 4°C to 14°C. Note the Desired Pressure in the table below depending on when the pressure is increased.
- Record the Perfusion Start Time (enzyme solution enters the pancreas) in the table below.
- Monitor temperature and pressure during pancreas perfusion and record in the table below.
- Optionally monitor the flow rate and record it in the table below.
- Stop perfusion after 10 minutes (12 minutes in the case of poor distension). If perfusion time exceeds 12 minutes, attach to this record a justification for the additional time.

Islets Lot Number: _____

Pancreas Perfusion Pressures & Temperatures

Continue to clean the pancreas during and after perfusion. Save all removed non-pancreatic tissue in the container from Section 5.11.

Post-perfusion trim finish time: _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 18 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

6.3 Final Trimmed Pancreas Weight

After perfusion and trimming are complete, weigh the additional tissue removed after the Initial Trimmed Pancreas Weight was determined (Section 5.7, above). Record this weight in row B of the table below, and calculate the Final Trimmed Pancreas Weight.

A. Initial Trimmed Pancreas Weight (from Section 5.7)	g
B. Additional Trimmed Tissue Weight	g
C. Final Trimmed Pancreas Weight ($A - B = C$)	g
D. Undigested Tissue Weight (from Section 7.3)	g
E. Digested Pancreas Tissue Weight ($C - D = E$)	g

Recorded by: _____ **Date:** _____

Verified by: _____ **Date:** _____

Determine the volume of CIT Enzyme Solution to be added to the Ricordi Digestion Chamber using the preparation table in the appropriate Attachment (B11, B13, B14) to SOP 3106.

Performed by: _____ **Date:** _____

- 6.4 Assemble the pancreas digestion equipment according to the institution's procedure. Use the 600 mL Ricordi Digestion Chamber (Biorep Technologies, Inc., Model No. 600-MUL-03 with screen WM-533, or Model No. 600-mDUR-03, with screen WM-533).

Performed by: _____ **Date:** _____

6.5 Pancreas Preparation for Digestion

Cut the pancreas into 5 to 15 similar sized pieces of 1 to 2.5 inches length and place the pieces in a Ricordi digestion chamber. Place 6 to 10 marbles into the digestion chamber and add CIT Enzyme Solution up to the point where the screen is to be placed. Place a 533 μ m woven stainless steel screen on top of the chamber and close it. Ensure that the digestion chamber is sealed properly to prevent leaking.

Performed by: _____ **Date:** _____

6.6 Pancreas Processing Times

Record information about the pancreas processing times in the table below. Calculate the Pancreas Preparation Time (Process Start Time, Section 5.2, to Perfusion Start Time, Section 6.2), and the Cold Ischemia Time (Cross Clamp Time, from donor records, to Perfusion Start Time, from Section 6.2) and record these in the table below.

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 19 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

	Date	Time
A. Cross Clamp (Donor Records)		
B. Process Start (Section 5.2)		
C. Perfusion Start (Section 6.2)		
	D. Pancreas Preparation Time (D = C – B)	Hours _____ Minutes _____
	E. Cold Ischemia Time* (E = C – A)	Hours _____ Minutes _____

*Cold Ischemia Time must be 12 hours or less. If the Cold Ischemia Time is more than 12 hours, immediately notify the site principal investigator.

Recorded by: _____ **Date:** _____

Calculate by: _____ **Date:** _____

Verified by: _____ **Date:** _____

If the site principal investigator is notified of excessive Cold Ischemia Time, complete the following:

Name of Person notified: _____

Notified by: _____

Date & Time Notified: _____, _____

7.0 ENZYMATIC PANCREAS DIGESTION

7.1 Pancreas Digestion

7.1.1 Add any remaining residual CIT Enzyme Solution to the recirculation flask for introduction into the digestion circuit.

Add 0 to 5 mL of Pulmozyme (2.5 mL/ampoule, 1 mg/mL) to the Ricordi Digestion Chamber

Volume of Pulmozyme (1 mg/mL) added: _____ mL

Performed by: _____ **Date:** _____

7.1.2 Start pumping the solution at a rate of 230 ± 20 mL/min to fill the system. Record this as the Digestion Start Time in the table in Section 7.2. Add as much CIT Digestion Solution to the recirculation flask as needed to fill the system and to completely eliminate air from the circuit.

Immediately begin recording the temperature inside the chamber, and the flow rate in the table in Section 7.2.

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 20 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Rock the chamber gently for the first 5 minutes and then decrease the flow rate to 110 ± 20 mL/min. Start shaking the chamber after 5 minutes. It takes approximately 3 - 5 minutes for the chamber to reach a target temperature of 32 to 38°C.

Verified by: _____ **Date:** _____

- 7.1.3 When tissue is observed in the circulating digest, take a 1 – 2 mL sample of the digest from the sampling port with a syringe. Place the digest sample in a 35 mm dish and add dithizone (DTZ) stain solution. Observe the digest under a microscope. Repeat this sampling (taking the same sample volume each time) and examination every 1-2 minutes during the digestion. Record the digestion chamber temperature, the flow rate and your observations on the stained sample in the table below. Maintain temperature between 32°C and 38°C, based on digest quality, considering the following factors that help in determining when to stop digestion and start dilution:

Factors	Ideal Ranges for Switching from Digestion to Dilution*
Amount of Tissue	3 to 6
Number of Islets	> 45 islets
% Free Islets	> 50%
% Fragmented (Over-digested) Islets	< 10%

*See definitions in Note, below.

Verified by: _____ **Date:** _____

Note:

Criteria for evaluating the digest and determining the end of digestion

- Estimate the amount of tissue by centering the tissue in the dish, viewing the mass with a microscope at 40X power, and estimating the amount of the visual field covered (6 = tissue covers entire visual field, 3 = tissue covers about 1/2 of the visual field, 0 = no tissue).
- Estimate the number of islets (a rough visual count, 10 – 20, 30 – 50, 80 – 90 islets, etc.).
- Estimate the % free islets (free islets versus the total number of islets, 25%, 50%, 90%, etc.). Free islets have less than 25% of the border attached to acinar tissue.
- Estimate the % fragmented islets (number of fragmented islets versus the total number of islets, 10%, 15%, 50%, etc.). Fragmented islets are those with a ragged border due to damage by overexposure to the enzyme (Over-digested).

- 7.1.4 When the decision to stop digestion is made, start dilution and collection of islets. Record the Dilution Start Time (= Digestion Stop Time) at the end of the table in Section 7.2 and calculate the Total Digestion Time.

Decided by: _____ **Date:** _____

Verified by: _____ **Date:** _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 21 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

7.2 Dilution and Collection of Islets

- Adjust the flow rate to 230 ± 20 mL/min, and continue shaking the digestion chamber.
- Add fresh RPMI 1640 at room temperature to the intake container as needed.
- Adjust the temperature of the chamber to ≤ 30 °C during dilution and collection.
 - If a large number of imbedded islets are observed in the digest, the chamber temperature may be maintained between 30°C and 38°C during dilution.
- Collect the digest into the 1L containers prepared in 3.3.2.
- Gently swirl each container periodically as it fills. When it reaches a volume of 1L, immediately decant the solution into 250 mL conical tubes for centrifugation at 170 X g and 2°C to 8°C for 3 to 4 minutes.
- Periodically take 1 to 2 mL samples of the diluted digest from the sample port with a syringe. Stain with Dithizone (DTZ) solution and observe the stained sample under a microscope. Record your observations in the table below.
- When no islets are observed in the stained samples and little tissue remains in the chamber, discontinue the addition of media to the system, collect the media remaining in the system, and stop the circulation pump.
- Record the Dilution Stop Time at the end of the table below, and calculate and record the Total Dilution Time.

Verified by: _____

Date: _____

Islets Lot Number: _____

Pancreas Digestion Record

Record Desired Temperatures and Desired Flow Rates in vacant cells based on Digestion Stop Time.

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 23 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

- 7.3 Remove the undigested pancreas material from the digestion chamber, weigh it, record the weight below, and in the table in Section 6.3. Calculate the weight of digested tissue in the table in Section 6.3.

Examine the undigested pancreas material remaining in the digestion chamber, and estimate the percentages of pancreatic tissue and connective tissue (should equal 100%). Record these estimates below.

Weight of undigested tissue remaining in chamber (record also in Section 6.3): _____ g
Calculate the Digested Pancreas Weight in Section 6.3 table, above.

Estimate of undigested pancreatic tissue: _____ %

Estimate of undigested connective tissue: _____ %

Performed by: _____ **Date:** _____

- 7.4 Tissue Recovery and Washing

- 7.4.1 Prior to the end of digestion prepare CIT Purification Solution and CIT Wash Solution according to DAIT SOP 3106, B02, and B12, respectively. Attach the record of preparation to this Production Batch Record and keep both solutions at 2°C to 8°C until used.
- 7.4.2 As tissue is collected during dilution, transfer it to 250 mL conical tubes for the first four liters and centrifuge at 170 X g and 2°C to 8°C for 3 to 4 minutes, to pellet the tissue.
- 7.4.3 Decant all of the supernatant and transfer pellets to a 1 L container containing 900 mL of CIT Wash Solution (keep cold).

NOTE: Be sure the flask is kept level during recombination to avoid tissue aggregation and hypoxic conditions.

- 7.4.4 If residual tissue remains, wash it with 3 to 5 mL of CIT Wash Solution.
- 7.4.5 After dilution is completed and all the tissue has been recombined into the CIT Wash Solution, mix the flask thoroughly by gentle swirling and transfer the contents into as many 250 mL sterile conical tubes as required. Centrifuge each tube at 170 X g and 2°C to 8°C for 3 to 4 minutes.
- 7.4.6 Wash the recombined tissue with CIT Wash Solution until the extracellular debris and DNA strings have been minimized. As the washing progresses, reduce the number of conical tubes to two, then one by combining tissue.

NOTE: If, during collection, DNA stings are observed after centrifugation with loose pellet formation, transfer the suspension portion of those tubes containing the majority of cells into one separate 250 mL conical tube, and keep it lying flat on the bench for 5 minutes after adding up to 200 mL of CIT Wash Solution and 200 µL (1 µg/mL) of Pulmozyme. After re-centrifugation, when the DNA strings have disappeared, recombine with other pellets.

- 7.4.7 After the washing is complete, centrifuge the final tube at 170 X g and 2°C to 8°C for 3 to 4 minutes and visually estimate the total packed tissue volume in the final 250 mL container. Aspirate the supernatant down to the pellet.

Total Packed Tissue Volume: _____ mL

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 24 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

- 7.4.8 Re-suspended the islets to 100 to 250 g or mL, depending on the amount of tissue, with CIT Purification Solution. Ensure there are no clumps (dissolve if necessary). Record the volume or weight.

Total Suspension Volume or Weight: _____ mL or _____ g

Verified by: _____ Date: _____

7.5 Pre-purification Islets Count

- 7.5.1 Re-suspend tissue evenly. Take two 100 μ L samples and count each sample once.

- 7.5.2 Perform pre-purification count according to the institution's procedure and record the data in the table below and attach spreadsheet, if used, to Production Batch Record.

Pre-purification Islets Counts & Calculations

Sample Volume	μ L			
Total Volume	mL			
Dilution Factor				
Diameter (μm), Factor	Counts		IPN (Avg.)	IEQ
50 – 100, 0.167				
101 – 150, 0.648				
151 – 200, 1.685				
201 – 250, 3.500				
251 – 300, 6.315				
301 – 350, 10.352				
> 350, 15.833				
		Sample Total		
		Suspension Total		
% Trapped				
% Fragmented				
Technicians' Initials				

Comments: _____

Verified by: _____ Date: _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 25 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

- 7.5.3 The maximum tissue volume for purification is 25 mL per COBE run. If the tissue volume is < 25 mL, centrifuge the islets suspension and re-suspend the tissue in 100 mL of CIT Purification Solution. If the tissue volume is > 25 mL, using the Packed Tissue Volume from Section 7.4.8, calculate the number of COBE runs required to process \leq 25 mL of packed tissue per run. Divide the tissue evenly into separate sterile 250 mL conical tubes and fill each to the 100 mL mark with additional CIT Purification Solution. During purification of the first tube, the additional conical tubes should be kept in the cold room or refrigerator for subsequent COBE runs (keep tube lying flat and mix occasionally to avoid tissue aggregation) until ready to be loaded into the COBE.

Number of conical tubes and COBE runs: _____

Volume of tissue distributed into each tube: _____ mL

Calculated by: _____ **Date:** _____

Verified by: _____ **Date:** _____

- 7.5.4 When ready to load the first COBE run, add 20 mL of Albumin Human USP, 25% Solution to the tissue and mix well. Continue to Section 8.2.11.

For subsequent COBE runs, centrifuge the conical tube at 170 X g and 2°C to 8°C for 3 – 4 minutes. Remove the supernatant, add 20 mL of Albumin Human USP, 25% Solution to the tissue and mix well to re-suspend. Bring the tissue suspension to 120 mL in a 250 mL tube or beaker with CIT Purification Solution. Continue to Section 8.2.11.

8.0 ISLETS PURIFICATION

8.1 COBE 2991 Preparation

Set up the COBE according to the Operational Manual and the institution's procedures. The COBE must be refrigerated or placed in a cold room.

- Prepare High (1.10 g/mL) and Low (1.06 g/mL) CIT Purification Density Gradients according to SOP 3106, B10, and file the records of their preparation with this Production Batch Record.
- Label 13 X 250 mL conical tubes with the COBE run number, and "W1" and fraction numbers 1 through 12 (See tables in Section 8.3). Label a 14th 250 mL conical tube with the COBE run number and "Bag."
- Fill tubes 1 through 12 with 225 mL of CMRL 1066, Supplemented, and store at 2°C to 8°C.

Verified by: _____ **Date:** _____

8.2 COBE 2991 Procedure – Gradient and Tissue Loading

- 8.2.1 Assemble the COBE bag onto COBE cell processor according to institution's procedure. Place clamps near the main line on all colored tubing except one line to be used for loading the COBE bag.

- 8.2.2 Place gradient-maker on magnetic stir plate and aseptically connect one end of size 16 tubing to gradient-maker and the other end to green tubing of the COBE bag.

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 26 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

- 8.2.3 Place a sterile stir bar into the left chamber (next to outlet) and turn on the stir plate.
- 8.2.4 Run tubing through pump and set pump to 60 mL/min.
- 8.2.5 Sanitize the exterior of all solution bottles before placing in the hood.
- 8.2.6 Pour 120 mL of the High Density Gradient (1.10 g/mL) into the left chamber of the gradient maker.
- 8.2.7 Start to pump High Density Gradient (1.10 g/mL) into COBE bag. Once this gradient reaches the bag, start the COBE at 1800 – 2000 rpm.
- 8.2.8 Once the entire 120 mL of High Density Gradient (1.10 g/mL) is loaded, remove excess air from the COBE bag by pressing Superout while unclamping the red tubing. Press the Hold button once the Bottom Gradient has reached the T (junction of red/green tube). Re-clamp the red tubing line and press the Stop/Reset button.
- 8.2.9 Wait for the final centrifugation of the digest tissue and then begin loading the continuous density gradient into the COBE bag (Section 7.5.4).
- Pour 125 mL High Density Gradient (1.10 g/mL) in the left chamber (nearest the outlet) of the gradient maker. Open and close the port between the two chambers just enough to fill the opening.
 - Pour 125 mL Low Density Gradient (1.06 g/mL) in the right chamber of gradient maker (away from outlet)
 - Start the COBE and ensure that the centrifuge speed is between 1800 and 2000 rpm.

Centrifuge Speed: _____ rpm

Recorded by: _____ Date: _____

- Open the port between the chambers, set pump to 20 mL/min and load gradient up to the T of the COBE bag tubing. Stop the pump when the gradient has reached the T-connection.

NOTE: Observe the gradient maker to ensure that gradients are mixing during the continuous gradient loading.

- 8.2.10 Load the continuous gradient by unclamping the green tubing and starting the pump. Load the entire 250 mL of continuous gradient at 20 mL/minute.
- 8.2.11 When all of the gradient has been loaded, stop the pump just as the last portion of the gradient enters the tubing attached to the gradient maker.

NOTE: COBE must remain spinning during the rest of the purification process. If abnormal signs appear from rotating seal (e.g. leak, unusual noise, burnt smell, etc.), replace COBE bag and make new density gradients.

- 8.2.12 Aseptically remove the tubing from gradient maker port and move it to the beaker with tissue. Reverse the pump to purge the air.
- 8.2.13 Load the tissue with the pump at a setting of 20 mL/min. Gently swirl the beaker to keep the tissue well-suspended during the loading.

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 27 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

- 8.2.14 To ensure tissue does not back-up on the gradient (a heavy tissue line observed on the gradient), periodically turn the pump off allowing tissue to enter the gradient and then turn the pump back on again. Repeat as necessary every 1 to 2 minutes.

NOTE: As an alternate, turn the pump off for 30 seconds, followed by loading tissue for 45 seconds.

- 8.2.15 As soon as the tissue is loaded, add 30 mL of additional CIT Purification Solution to the 250 mL beaker to rinse. Load this rinse onto the COBE.
- 8.2.16 After the last portion of the rinse has entered the COBE bag, stop the pump.
- 8.2.17 Vent the system by carefully unclamping the red tubing. Re-clamp the tubing when liquid (capping solution) is approximately one inch above the ceramic seal. This is the start of centrifugation time.

NOTE: Air left in the ceramic rotating seal can cause seal failure which may lead to leaking, seal occlusion and possible system shutdown due to overpressure during Superout.

- 8.2.18 Clamp the green line and allow the COBE to spin for 3 minutes. Record data on Purification Data Log for each COBE run, below.

Verified by: _____ **Date:** _____

8.3 COBE 2991 Procedure – Tissue Collection

- 8.3.1 During the 3 minute spin disconnect tubing from the pump. Prepare for collection of tissue fractions.
- 8.3.2 Verify that the Superout Rate is set at 100 mL/min.
- 8.3.3 After 3 minute spin slowly remove the blue clamp on the green line and quickly press the Superout button.
- 8.3.4 Collect the first 150 mL of effluent into the conical tube labeled “W” and 12 X 25 mL fractions into the numbered conical tubes each pre-filled with 225 mL CMRL 1066, Supplemented, CIT Modifications, as described on the Purification Data Log for each respective COBE run.
- 8.3.5 Once the fractions are collected, stop the COBE and aseptically collect the contents of the COBE bag into a 250 mL conical tube labeled “bag.” Discard the COBE bag and tubing.
- 8.3.6 Dilute the COBE bag contents up to 200 mL with CMRL 1066, Supplemented, CIT Modifications. Take a 200 µL sample and place it into 35 mm dish. Stain the sample with dithizone according to the institution’s procedure and examine it for the presence of islets. If a significant number of free islets are present keep the diluted COBE bag contents at 2°C to 8°C for further processing as instructed in Section 8.4.1. If there are not a significant number of free islets, discard the COBE bag contents.
- 8.3.7 To evaluate each COBE fraction quickly, gently but thoroughly mix each fraction from Section 8.3.4, then quickly transfer a 0.5 mL sample to one well of a 12-well microtiter plate and 0.5 mL of the W fraction to a 35 mm dish.
- 8.3.8 Stain each sample with dithizone according to the institution’s procedure and observe for islets. Record Islets Purity (%) and disposition of each fraction on the Purification Data Log for each COBE run.

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 28 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

8.3.9 Centrifuge the 250 mL tubes for 3 minutes at 140 X g and 2°C to 8°C. Record Packed Tissue Volumes of each COBE fraction on the Purification Data Log for each respective COBE run. Discard supernatant.

8.3.10 Combine the islets fractions by transferring the pellets with 10 mL pipets into four labeled 250 mL conical tubes containing 100 mL of CMRL 1066, Supplemented, to obtain the following purity levels after recombination:

- High Purity ($\geq 70\%$) (H),
- Middle Purity (40% to 69%) (M),
- Low Purity (30% to 39%) (L), and
- Supplementary Purification Islets ($< 30\%$) (S).

Discard fractions (D) that contain little or no tissue. For the other four categories of islets purity, keep the conical tubes flat on the bench at room temperature until the tissue of all COBE runs has been combined into the respective conical tubes.

NOTE: Depending on the analysis and disposition of each fraction, there may be up to one 250 mL conical tube for each Purity Level (High, Middle, Low Purity Islets), and one 250 mL conical tube for the Supplementary Purification Islets, if there are any.

8.3.11 Repeat steps 8.2.1 to 8.3.10 for each COBE purification run. Combine fractions of similar purity into the 250 mL conical tubes prepared in Section 8.3.10.

NOTE: **Scoring Guidelines for purified layers in Purification Data Logs:**

- Packed Tissue Volume: estimate of the tissue volume in the individual conical tubes after they have centrifuged for 3 minutes at 140 X g and 2°C to 8°C.
- % Purity: estimate relative amount (%) of islets to total tissue.
- H M L S D: This is the disposition of each fraction as defined in the column header.

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 29 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Repeat this purification process for each of the tubes.

Purification Data Log, COBE Run #1:

Layer	Medium	Amount
Capping Layer	CIT Purification Solution	30 mL
Tissue Layer	_____ mL packed tissue in this COBE Run, plus 20 mL of Albumin Human USP, 25% Solution, and q.s. to 120 g with CIT Purification Solution	120 g
Density Gradients	Low Density Gradient (1.06 g/mL)	125 g
	High Density Gradient (1.10 g/mL)	125 g
Bottom	High Density Gradient (1.10 g/mL)	120 g
Centrifuge Start Time		Centrifuge Stop Time

#	CMRL 1066, Supplemented Pre-fill Vol. (mL)	Fraction Volume Collected (mL)	Packed Tissue Volume (mL)	Comments	Islet Purity (%)	Disposition: H: High, M: Middle, L: Low, S: Supplementary, D: Discard (Circle One)
W	0	150 mL				H M L S D
1	225	25				H M L S D
2	225	25				H M L S D
3	225	25				H M L S D
4	225	25				H M L S D
5	225	25				H M L S D
6	225	25				H M L S D
7	225	25				H M L S D
8	225	25				H M L S D
9	225	25				H M L S D
10	225	25				H M L S D
11	225	25				H M L S D
12	225	25				H M L S D
Bag	0	95				S D

Comments on purification: _____

Recorded by: _____

Date: _____

Verified by: _____

Date: _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 30 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Purification Data Log, COBE Run #2

Layer	Medium	Amount
Capping Layer	CIT Purification Solution	30 mL
Tissue Layer	_____ mL packed tissue in this COBE Run, plus 20 mL of Albumin Human USP, 25% Solution, and q.s. to 120 g with CIT Purification Solution	120 g
Density	Low Density Gradient (1.06 g/mL)	125 g
Gradients	High Density Gradient (1.10 g/mL)	125 g
Bottom	High Density Gradient (1.10 g/mL)	120 g
Centrifuge Start Time		Centrifuge Stop Time

#	CMRL 1066, Supplemented Pre-fill Vol. (mL)	Fraction Volume Collected (mL)	Packed Tissue Volume (mL)	Comments	Islet Purity (%)	Disposition: H: High, M: Middle, L: Low, S: Supplementary, D: Discard (Circle One)
W	0	150				H M L S D
1	225	25				H M L S D
2	225	25				H M L S D
3	225	25				H M L S D
4	225	25				H M L S D
5	225	25				H M L S D
6	225	25				H M L S D
7	225	25				H M L S D
8	225	25				H M L S D
9	225	25				H M L S D
10	225	25				H M L S D
11	225	25				H M L S D
12	225	25				H M L S D
Bag	0	95				S D

Comments on purification: _____

Recorded by: _____

Date: _____

Verified by: _____

Date: _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 31 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Purification Data Log, COBE Run #3

Layer	Medium	Amount
Capping Layer	CIT Purification Solution	30 mL
Tissue Layer	_____ mL packed tissue in this COBE Run, plus 20 mL of Albumin Human USP, 25% Solution, and q.s. to 120 g with CIT Purification Solution	120 g
Density Gradients	Low Density Gradient (1.06 g/mL)	125 g
	High Density Gradient (1.10 g/mL)	125 g
Bottom	High Density Gradient (1.10 g/mL)	120 g
Centrifuge Start Time		Centrifuge Stop Time

#	CMRL 1066, Supplemented Pre-fill Vol. (mL)	Fraction Volume Collected (mL)	Packed Tissue Volume (mL)	Comments	Islet Purity (%)	Disposition: H: High, M: Middle, L: Low, S: Supplementary, D: Discard (Circle One)
W	0	150				H M L S D
1	225	25				H M L S D
2	225	25				H M L S D
3	225	25				H M L S D
4	225	25				H M L S D
5	225	25				H M L S D
6	225	25				H M L S D
7	225	25				H M L S D
8	225	25				H M L S D
9	225	25				H M L S D
10	225	25				H M L S D
11	225	25				H M L S D
12	225	25				H M L S D
Bag	0	95				S D

Comments on purification: _____

Recorded by: _____

Date: _____

Verified by: _____

Date: _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 32 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Purification Data Log, COBE Run #4

Layer	Medium	Amount
Capping Layer	CIT Purification Solution	30 mL
Tissue Layer	_____ mL packed tissue in this COBE Run, plus 20 mL of Albumin Human USP, 25% Solution, and q.s. to 120 g with CIT Purification Solution	120 g
Density	Low Density Gradient (1.06 g/mL)	125 g
Gradients	High Density Gradient (1.10 g/mL)	125 g
Bottom	High Density Gradient (1.10 g/mL)	120 g
Centrifuge Start Time		Centrifuge Stop Time

#	CMRL 1066, Supplemented Pre-fill Vol. (mL)	Fraction Volume Collected (mL)	Packed Tissue Volume (mL)	Comments	Islet Purity (%)	Disposition: H: High, M: Middle, L: Low, S: Supplementary, D: Discard (Circle One)
W	0	150				H M L S D
1	225	25				H M L S D
2	225	25				H M L S D
3	225	25				H M L S D
4	225	25				H M L S D
5	225	25				H M L S D
6	225	25				H M L S D
7	225	25				H M L S D
8	225	25				H M L S D
9	225	25				H M L S D
10	225	25				H M L S D
11	225	25				H M L S D
12	225	25				H M L S D
Bag	0	95				S D

Comments on purification: _____

Recorded by: _____

Date: _____

Verified by: _____

Date: _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 33 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Purification Data Log, COBE Run #5

Layer	Medium	Amount
Capping Layer	CIT Purification Solution	30 mL
Tissue Layer	_____ mL packed tissue in this COBE Run, plus 20 mL of Albumin Human USP, 25% Solution, and q.s. to 120 g with CIT Purification Solution	120 g
Density Gradients	Low Density Gradient (1.06 g/mL)	125 g
	High Density Gradient (1.10 g/mL)	125 g
Bottom	High Density Gradient (1.10 g/mL)	120 g
Centrifuge Start Time		Centrifuge Stop Time

#	CMRL 1066, Supplemented Pre-fill Vol. (mL)	Fraction Volume Collected (mL)	Packed Tissue Volume (mL)	Comments	Islet Purity (%)	Disposition: H: High, M: Middle, L: Low, S: Supplementary, D: Discard (Circle One)
W	0	150				H M L S D
1	225	25				H M L S D
2	225	25				H M L S D
3	225	25				H M L S D
4	225	25				H M L S D
5	225	25				H M L S D
6	225	25				H M L S D
7	225	25				H M L S D
8	225	25				H M L S D
9	225	25				H M L S D
10	225	25				H M L S D
11	225	25				H M L S D
12	225	25				H M L S D
Bag	0	95				S D

Comments on purification: _____

Recorded by: _____

Date: _____

Verified by: _____

Date: _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 34 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Note: If the initial purification process, above, did not yield a sufficient number of sufficiently pure islets for transplant, and there is a substantial quantity of tissue containing impure islets in the Middle and/or Low Purity Islets 250 mL conical tubes, and/or in the Supplementary Purification 250 mL conical tube, follow the procedure in Section 8.4, below.

8.4 Supplementary Purification Fractions and COBE Bag Contents Processing

8.4.1 If, upon examination of the COBE bag contents, a significant number of islets is present (See Section 8.3.6), centrifuge the 250 mL conical tube containing the diluted COBE bag contents at 140 X gravity and 2°C to 8°C for three minutes, and transfer the packed tissue to the Supplementary Purification Islets 250 mL conical tube.

8.4.2 List all fractions combined for Supplementary Purification:

COBE Run #	Fractions and/or COBE Bags Combined for Supplementary Purification
1	
2	
3	
4	
5	

Recorded by: _____ **Date:** _____

Verified by: _____ **Date:** _____

8.4.3 Bring the volume of the Supplementary Purification Islets 250 mL conical tube to 100 to 250 mL with CMRL 1066, Supplemented, CIT Modifications, and take one or two 100 µL samples for counting, if desired.

8.4.4 Dilute the Supplementary Purification Islets to 250 mL with CMRL 1066, Supplemented, CIT Modifications. Lay the tube on its side at 2°C to 8°C if counts are performed.

Verified by: _____ **Date:** _____

8.4.5 If desired, count islets according to the institution's procedure in the Supplementary Purification Islets sample and record counts in the table below and attach any spreadsheets used. Indicate in the Comments space if the tissue will be re-purified. Supplementary Purification may be indicated if there are a significant number of islets (greater than 50,000 IEQ). If Supplementary Purification is to be performed, record which of the two procedures will be used on the Comments lines below the Counts table, and proceed to Section 9.0. If Supplementary Purification is not to be performed, record the disposition of the Supplementary Purification Islets on the Comments lines below the Counts table.

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 35 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Optional Pre-supplementary Purification Islets Counts & Calculations

Sample Volume	μL			
Total Volume	mL			
Dilution Factor				
Diameter, Factor	Counts		IPN (Avg.)	IEQ
50 – 100, 0.167				
101 – 150, 0.648				
151 – 200, 1.685				
201 – 250, 3.500				
251 – 300, 6.315				
301 – 350, 10.352				
> 350, 15.833				
		Sample Total		
		Suspension Total		
% Trapped				
% Fragmented				
Technicians' Initials				

Comments: _____

Recorded by: _____ Date: _____

Verified by: _____ Date: _____

Decided by: _____ Date: _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 36 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

8.5 Tissue Preparation for Re-purification

If the decision in Section 8.4, is to perform a Supplementary Purification of the islets, centrifuge the 250 mL conical tube containing all the Supplementary Purification Islets at 140 X gravity and 2°C to 8°C for three minutes. Remove and discard the supernatant.

Performed by: _____ **Date:** _____

Verified by: _____ **Date:** _____

9.0 ISLETS SUPPLEMENTARY PURIFICATION

If islets tissue insufficiently purified by the procedure described in Section 8.0 is present, this tissue may be re-purified by one of the three procedures defined in SOP 3109. Cross out all three references, if no Supplementary Purification is performed. Cross out the two references not used, if Supplementary Purification is performed.

9.1 SOP 3109, B01, Supplementary Purification, OptiPrep Procedure & Record

9.2 SOP 3109, B02, Supplementary Purification, Continuous Biocoll Procedure & Record

9.3 SOP 3109, B03, Supplementary Purification, Discontinuous Polysucrose Procedure & Record

File the Supplementary Purification record with this Production Batch Record.

Recorded by: _____ **Date:** _____

Approved by: _____ **Date:** _____

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 37 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

10.0 POST-PURIFICATION ISLETS COUNT

- 10.1 After all islets are combined into the three Purity Levels, wash each Purity Level once with CIT Culture Media prepared according to DAIT SOP 3106, B04. Allow the tissue in the conical tubes to settle for 3 to 5 minutes. After the tissue in each purity level has settled, remove the supernatant and re-suspend the final tissue in 50 to 250 mL of CIT Culture Media in T-75 flasks labeled for each Purity Level with Lot Number and isolation date.

Verified by: _____ Date: _____

- 10.2 Gently mix each Purity Level and take two 100 μ L samples of each for Post-purification Islet Count. Enter the count data in the table below, attach a spreadsheet, if used, and calculate the Total Islet Number (IPN) and Total IEQ. The contents of these T-75 flasks are now ready to proceed to Islet Culture, Section 11.

Sampled by: _____ Date: _____

Post-purification Islets Counts

	High Purity				Middle Purity				Low Purity			
Sample Volume	μ L				μ L				μ L			
Total Volume	mL				mL				mL			
Dilution Factor												
Diameter, Factor	Counts	Avg.	IEQ		Counts	Avg.	IEQ		Counts	Avg.	IEQ	
50 – 100, 0.167												
101 – 150, 0.648												
151 – 200, 1.685												
201 – 250, 3.500												
251 – 300, 6.315												
301 – 350, 10.352												
> 350, 15.833												
Total												
% Trapped												
% Fragmented												
% Purity												
Islet Quality Grade*												
Technicians' Initials												

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 38 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Post-purification Islets Calculations

	High Purity	Middle Purity	Low Purity	Total
Post-purification IPN				
Post Purification IEQ				
Pre-purification IEQ (Section 7.5.2)				
IEQ Recovery (%) (from Pre-purification IEQ)				
Total IEQ/g of Final Trimmed Pancreas (Section 6.3)				
Comments				

*See Note, below, for Islets Quality Grade guidelines

Calculated by: _____ Date: _____

Verified by: _____ Date: _____

Note: Islets Quality Grade

Grade the quality of the islets based on these parameters and criteria:

Parameter	0 Points	1 Point	2 Points
Shape (3D)	flat/planar	in between	spherical
Border (2D)	irregular	in between	well-rounded
Integrity	fragmented	in between	solid/compact
Single Cells	many	a few	almost none
Diameter	all < 100 µm	a few > 200 µm	> 10% > 200 µm

Add up the points for each sample to obtain the following grades:

- 9 to 10 points = A
- 7 to 8 points = B
- 4 to 6 points = C
- 2 to 3 points = D
- 0 to 1 point = F

Islets Lot Number: _____

Document No. SOP 3101, B02-1	Revision No. 07	Effective Date 06 Aug 2011	Supersedes Date 02 May 2011	Page 39 of 39
Document Title: PHPI, MASTER PRODUCTION BATCH RECORD (PRODUCT CODES PHPI-A-01 & PHPI-L-01)				

Will the CIT Culture Media, the CIT Transplant Wash Media, and the CIT Transplant Media for this lot contain no drug, or Lisofylline?

Circle one of the following:

Islets Alone (Clinical Protocols 03 – 07)

Islets with Lisofylline (Clinical Protocol 02 only)

Approved by: _____ Date: _____

If “Islets Alone” is circled above, continue recording the manufacturing process in Part 2A (SOP 3101, B02-2A).

If “Islets with Lisofylline” is circled above, continue recording the manufacturing process in Part 2B (SOP 3101, B02-2B).

Note: Part 1 of the Production Batch Record must be combined with either Part 2A or Part 2B for review and approval.

Islets Lot Number: _____