

**The NIH CIT Consortium Chemistry Manufacturing Controls Monitoring Committee:**

J. Ansite, A.N. Balamurugan, B. Barbaro, J. Battle, D. Brandhorst, J. Cano, X. Chen, S. Deng, D. Feddersen, A. Friberg, T. Gilmore, J.S. Goldstein, E. Holbrook, A. Khan, T. Kin, J. Lei, E. Linetsky, C. Liu, X. Luo, K. McElvaney, Z. Min, J. Moreno, D. O’Gorman, K.K. Papas, G. Putz, C. Ricordi, G. Szot, T. Templeton, L. Wang, J.J. Wilhelm, J. Willits, T. Wilson, X. Zhang

**The NIH CIT Consortium**

**Emory University:** J. Avila, B. Begley, J. Cano, S. Carpentier, E. Holbrook, J. Hutchinson, C.P. Larsen, J. Moreno, M. Sears, N.A. Turgeon, D. Webster

**Massachusetts General Hospital:** S. Deng, J. Lei, J.F. Markmann

**NIAID:** N.D. Bridges, C.W. Czarniecki, J.S. Goldstein, G. Putz, T. Templeton, T. Wilson

**NIDDK:** T.L. Eggerman

**Northwestern University:** P. Al-saden, J. Battle, X. Chen, A. Hecyk, H. Kissler, X. Luo, M. Molitch, N. Monson, E. Stuart, A. Wallia, L. Wang, S. Wang, X. Zhang

**University of Alberta, Edmonton:** D. Bigam, P. Campbell, P. Dinyari, T. Kin, N. Kneteman, J. Lyon, A. Malcolm, D. O’Gorman, C. Onderka, R. Owen, R. Pawlick, B. Richer, S. Rosichuk, D. Sarman, A. Schroeder, P.A. Senior, A.M.J. Shapiro, L. Toth, V. Toth, W. Zhai

**University of California–San Francisco:** K. Johnson, J. McElroy, A.M. Posselt, M. Ramos, T. Rojas, P.G. Stock, G. Szot

**University of Illinois, Chicago:** B. Barbaro, J. Martellotto, J. Oberholzer, M. Qi, Y. Wang

**University of Iowa (Data Coordinating Center):** L. Bayman, K. Chaloner, W. Clarke, J.S. Dillon, C. Diltz, G.C. Doelle, D. Ecklund, D. Feddersen, E. Foster, L. G. Hunsicker, C. Jaspersen, D-E Lafontant, K. McElvaney, T. Neill-Hudson, D. Nollen, J. Qidwai, H. Riss, T. Schwieger, J. Willits, J. Yankey

**University of Miami:** R. Alejandro, A.C. Corrales, R. Faradji, T. Froud, A.A. Garcia, E. Herrada, H. Ichii, L. Inverardi, N. Kenyon, A. Khan, E. Linetsky, J. Montelongo, E. Peixoto, K. Peterson, C. Ricordi, J. Szust, X. Wang

**University of Minnesota:** M.H. Abdulla, J. Ansite, A.N. Balamurugan, M.D. Bellin, M. Brandenburg, T. Gilmore, J. V. Harmon, B.J. Hering, R. Kandaswamy, G. Loganathan, K. Mueller, K.K. Papas, J. Pedersen, J.J. Wilhelm, J. Witson

**University of Pennsylvania:** C. Dalton-Bakes, H. Fu, M. Kamoun, J. Kearns, Y. Li, C. Liu, E. Luning-Prak, Y. Luo, E. Markmann, Z. Min, A. Naji, M. Palanjan, M. Rickels, R. Shlansky-Goldberg, K. Vivek, A.S. Ziaie

**University of Wisconsin:** L. Fernandez, D.B. Kaufman, L. Zitur

**Uppsala University:** D. Brandhorst, A. Friberg, O. Korsgren

*Supported by grants from the National Institute of Allergy and Infectious Diseases and the National Institute for Diabetes and Digestive and Kidney Diseases.*

- At Emory University, U01AI089317.
- At Northwestern University, U01AI089316.
- At the University of Alberta, Edmonton: U01AI065191.
- At the University of California, San Francisco, U01DK085531.
- At the University of Illinois, Chicago, 5U01DK070431-10.
- At the University of Iowa, U01DK070431.
- At the University of Miami, U01DK070460.
- At the University of Minnesota, U01AI065193.
- At the University of Pennsylvania, U01DK070430.
- At Uppsala University, U01AI065192.

*In addition, the study was supported by the following GCRC and CTSA awards:*

- At Emory University: UL1TR000454.
- At Northwestern University: 5UL1RR025741 and 8UL1TR000150.
- At the University of California, San Francisco, UL1TR000004.
- At the University of Illinois, Chicago, UL1TR000050.
- At the University of Miami: 1UL1TR000460.
- At the University of Minnesota: 5M01-RR000400 and UL1TR000114.
- At the University of Pennsylvania: UL1TR000003.

Address correspondence to: Camillo Ricordi MD, Chairman, CIT Steering Committee,  
ricordi@miami.edu

## To cite this article

*PHPI, MPBR, Part 2A (Product Code PHPI-A-01, Islets Alone)– Standard Operating Procedure of the NIH Clinical Islet Transplantation Consortium*

CellR4 2017; 5 (2): e2288

|  |                    |                                     |                                 |              |
|--|--------------------|-------------------------------------|---------------------------------|--------------|
| Document No.<br>SOP 3101, B02-2A   | Revision No.<br>04 | Effective Date<br>04 September 2009 | Supersedes Date<br>21 July 2009 | Page 1 of 32 |
| Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE) |                    |                                     |                                 |              |

**Note:** Use this document only if “Islets Alone” are being manufactured.

## 11.0 ISLET CULTURE

- 11.1 For product characterization tests samples, gently re-suspend the contents of the High Purity ( $\geq 70\%$ ) Islets culture flask. Based on the count results in Section 10, take a sample containing  $\geq 400$  IEQ for a Pre-culture Glucose Stimulated Insulin Release Test according to the institution's procedure. This islets sample is cultured in a culture dish simultaneously with, but separately from, the bulk islets product. Report Result in Section 14.5 and on the Interim Certificate of Analysis.

Also, take samples of the High Purity Islets suspension for the Pre-culture DNA Content, and Nuclei Measurement product characterization tests according to the table, below. Report the results of these tests in Section 17.3.

| CHARACTERIZATION TEST  | IEQ     | IEQ/mL | SAMPLE REMOVED (mL) |
|--|---------|--------|---------------------|
| Example –Low Yield   | 400     | 1,000  | 0.40 mL             |
| Example – High Yield   | 400     | 5,000  | 0.08 mL             |
| <b>Certificates of Analysis</b>                                |         |        |                     |
| REQUIRED PRE-CULTURE GLUCOSE STIMULATED INSULIN RELEASE        | 400     |        |                     |
| <b>Optional Product Characterization, For Information Only</b> |         |        |                     |
| PRE-CULTURE DNA CONTENT  | 3 X 100 |        |                     |
| PRE-CULTURE NUCLEI MEASUREMENT                                 | 3 X 100 |        |                     |
| Sampled by:  |         |        | Date:               |
| Verified by:   |         |        | Date:               |

- 11.2 Calculate the number of T-175 culture flasks needed for a target of 20,000 to 30,000 IEQ/Flask using the equation (Round decimals up to the next higher whole number of flasks):

$$\frac{\text{IEQ in Purity Level}}{(\text{20,000 to 30,000 IEQ/Flask}) \times \text{Purity (in decimal form)}} = \# \text{ of T-175 Culture Flasks}$$

| Purity Level            | IEQ Counts | Purity | Target IEQ/Flask | Number of T-175 Culture Flasks |
|-------------------------|------------|--------|------------------|--------------------------------|
| Example – High Purity   | 352,423    | 0.95   | 27,500           | 13.48988, rounded up to 14     |
| Example – Middle Purity | 53,817     | 0.50   | 25,000           | 4.30536 rounded up to 5        |
| High Purity             |            |        |                  |                                |
| Middle Purity           |            |        |                  |                                |
| Low Purity              |            |        |                  |                                |
| Calculated by:          |            |        | Date:            |                                |
| Verified by:            |            |        | Date:            |                                |

|   |                           |  |  |                     |
|---|---------------------------|--|--|---------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 2 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                     |

- 11.3 Obtain the calculated number of sterile T-175 flasks, inspect each for cracks, and label them.

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

- 11.4 Transfer 20,000 to 30,000 IEQ to each T-175 culture flask and bring the volume to 30 mL with CIT Culture Media

| <b>Fraction</b>           | <b>Number of T-175 Culture Flasks</b> | <b>Media Needed (10 mL/flask)</b> | <b>CIT Culture Media Volume (Section 10.2)</b> | <b>CIT Culture Media Added or Removed</b> |
|---------------------------|---------------------------------------|-----------------------------------|--|---|
| Example 1 – High Purity   | 14                                    | 140 mL                            | 100 mL   | + 40 mL                                   |
| Example 2 – Middle Purity | 5                                     | 150 mL                            | 120 mL   | + 30 mL                                   |
| Example 3 – Low Purity    | 2                                     | 60 mL                             | 100 mL   | – 40 mL                                   |
| High Purity               |                                       |                                   |  |   |
| Middle Purity             |                                       |                                   |  |   |
| Low Purity                |                                       |                                   |  |   |
| Calculated by:            |                                       |                                   | Date:  |   |
| Verified by:              |                                       |                                   | Date:  |   |
| Performed by:             |                                       |                                   | Date:  |   |

- 11.5 Add 15 mL of CIT Culture Media to the culture dish containing the sample for Glucose Stimulated Insulin Release Assay (Section 11.1) and culture its contents with the High Purity Islets.

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

- 11.6 Place all the flasks of High Purity Islets in an incubator at 37°C, 95% air, and 5% carbon dioxide and record the date and time.

Date and time High Purity Islets flasks placed in 37°C incubator: \_\_\_\_\_

Record this date and time in the table in Section 12.5.

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

The islets' culture must end (Section 12.5) between 36 and 72 hours of the start time. Calculate these dates and times.

Date and time of minimum culture: \_\_\_\_\_

Date and time of maximum culture: \_\_\_\_\_

**Calculated by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

|   |                           |  |  |                     |
|---|---------------------------|--|--|---------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 3 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                     |

Notify the Site Principal Investigator, or designee, of these dates and times.

**Name of Person notified:** \_\_\_\_\_

**Notified by:** \_\_\_\_\_ **Date & Time Notified:** \_\_\_\_\_

- 11.7 Place all the flasks of Middle and Low Purity Islets in an incubator at 22°C, 95% air, and 5% carbon dioxide with the T-neck in the up position and record the date and time.

Date and time Middle and Low Purity Islets flasks placed in 22°C incubator: \_\_\_\_\_

Record this date and time in the table in Section 12.5.

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

- 11.8 Media Change

- 11.8.1 After 12 to 24 hours remove all the flasks from the incubators and record the date(s) and time(s) that each purity level of islet product is removed from the incubator(s) in the table in Section 12.5.

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

- 11.8.2 Inspect the contents of each flask for gross appearance, cloudiness, stranding or clumping. Using a microscope, examine the morphology of the islets, including the extent of fragmentation and the numbers of single cells; and the fluid in each flask for microorganisms. Signs of contamination (cloudiness, microorganisms upon microscopic examination) or unusual islets morphology, including extensive fragmentation or large numbers of single cells, must be reported to the Site Principal Investigator, or designee, immediately, and investigated according to the institution's procedures. Record observations and dispositions of flasks below.

---



---



---



---



---



---



---

**Inspected by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

If the Site Principal Investigator, or designee, is notified of any unusual islets morphology or evidence of microbial contamination, complete the following:

**Name of Person notified:** \_\_\_\_\_

**Notified by:** \_\_\_\_\_

**Date & Time Notified:** \_\_\_\_\_



|   |                           |  |  |                     |
|---|---------------------------|--|--|---------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 4 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                     |

- 11.8.3 Equilibrate the CIT Culture Media at room temperature. Place each flask in the BSC, tilt each at a 45° angle, and allow the islets to settle for 2 to 3 minutes. Aseptically remove 20 mL of supernatant media from each flask, and place all the removed supernatant from each purity level in as many containers as necessary for that purity level.

Add 20 mL of fresh CIT Culture Media to each flask, and replace the cap on each flask.

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

- 11.8.4 Transfer the supernatants to 250 mL conical tubes and centrifuge at 140 X g for 3 minutes. Remove supernatant and transfer tissue (if present) to a separate T-175 culture flask for each purity level.

|                                       | <b>High Purity Supernatant</b> |           | <b>Middle Purity Supernatant</b> |           | <b>Low Purity Supernatant</b> |           |
|---------------------------------------|--------------------------------|-----------|----------------------------------|-----------|-------------------------------|-----------|
| <b>Tissue Observed and recovered?</b> | <b>Yes</b>                     | <b>No</b> | <b>Yes</b>                       | <b>No</b> | <b>Yes</b>                    | <b>No</b> |

**Checked by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

If no tissue is observed, discard the supernatant as biohazardous waste.

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

- 11.9 Place all the T-175 culture flasks (High, Middle, and Low Purity Levels) into an incubator at 22°C, 95% air, and 5% carbon dioxide with the T-neck in the up position, and record the date(s) and time(s) that each purity level of islet product is placed in the incubator(s) in the table in Section 12.5.

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

## 12.0 ISLET PREPARATION FOR TRANSPLANT

- 12.1 Record the date and time scheduled for transplant of this lot of islets.

Scheduled Islet Transplant Date: \_\_\_\_\_

Scheduled Islet Transplant Time: \_\_\_\_\_

**Recorded by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

- 12.2 Physician's Order for Transplant

Verify that the physician's signed order for transplant (if used by the institution) is present, and the order, or a copy, is attached to this batch record.

Yes No (Circle One)

Physician's Name: \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

|   |                           |  |  |                     |
|---|---------------------------|--|--|---------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 5 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                     |

### 12.3 Recipient & Donor Information

From the source documents record the information about the prospective recipient in the table below. Attach a copy of the Request for Islet Transplant form to this Production Batch Record.

|                                     | <b>Islet Recipient Information</b> | <b>Donor Information</b> |
|-------------------------------------|------------------------------------|--------------------------|
| Hospital Name                       |                                    | UNOS or DDD #            |
| Recipient Medical Record Number     |                                    |                          |
| Recipient Study ID #                |                                    |                          |
| Date of Birth                       |                                    |                          |
| Gender                              |                                    |                          |
| ABO                                 |                                    |                          |
| CMV Status                          |                                    |                          |
| Allergies (Cipro, Penicillin, etc.) |                                    |                          |
| Current Weight (kg)                 |                                    |                          |

**Recorded by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

Compare this information with the Donor information in Section 4.4.

Blood Type Compatible? Yes No (Circle One)

CMV Status Reviewed? Yes No (Circle One)

Allergies Reviewed? Yes No (Circle One)

Information Reviewed with Clinician? Yes No (Circle One)

**Compared by:** \_\_\_\_\_ **Date:** \_\_\_\_\_  
Lab Manager or designee

**Reviewed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

### 12.4 Before the scheduled transplant time:

12.4.1 Prepare the laboratory and the Biological Safety Cabinets (BSCs) according to the institution's procedure.

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

12.4.2 In a BSC prepare CIT Transplant Wash Media and CIT Transplant Media according to DAIT SOP 3106, B05 and B06, respectively, and attach the record of preparation to this Production Batch Record. Equilibrate these media to room temperature before use.

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

|   |                           |  |  |                     |
|---|---------------------------|--|--|---------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 6 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                     |

12.5 Remove all the islet product flasks from the incubator(s) and record the date and time in the table below.

|  |      | High Purity<br>Flasks | Middle Purity<br>Flasks | Low Purity<br>Flasks | Recorded<br>by | Verified<br>by |
|--|------|-----------------------|-------------------------|----------------------|----------------|----------------|
| <b>1<sup>st</sup> Culture Start</b>        | Date |                       |                         |                      |                |                |
|  | Time |                       |                         |                      |                |                |
| <b>1<sup>st</sup> Culture Stop</b>         | Date |                       |                         |                      |                |                |
|  | Time |                       |                         |                      |                |                |
| <b>1<sup>st</sup> Culture Time (Hours)</b> |      |                       |                         |                      |                |                |
| <b>2<sup>nd</sup> Culture Start</b>        | Date |                       |                         |                      |                |                |
|  | Time |                       |                         |                      |                |                |
| <b>2<sup>nd</sup> Culture Stop</b>         | Date |                       |                         |                      |                |                |
|  | Time |                       |                         |                      |                |                |
| <b>2<sup>nd</sup> Culture Time (Hours)</b> |      |                       |                         |                      |                |                |
| <b>Total Culture Time (Hours)</b>          |      |                       |                         |                      |                |                |

Is the date and time of the 2<sup>nd</sup> Culture Stop within the dates and times of minimum and maximum culture calculated in Section 11.6?

Yes                      No                      (Circle One)

If it is not, immediately notify the Principal Investigator, or designee.

**Recorded by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

If the Site Principal Investigator, or designee, is notified, complete the following:

**Name of Person notified:** \_\_\_\_\_

**Notified by:** \_\_\_\_\_ **Date & Time Notified:** \_\_\_\_\_



|   |                           |  |  |                     |
|---|---------------------------|--|--|---------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 7 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                     |

- 12.6 Inspect the contents of each flask for gross appearance, cloudiness, stranding or clumping. Using a microscope, examine the morphology of the islets, including the extent of fragmentation and the numbers of single cells; and the fluid in each flask for microorganisms. Signs of contamination (cloudiness, microorganisms upon microscopic examination) or unusual islets morphology, including extensive fragmentation or large numbers of single cells, must be reported to the Site Principal Investigator, or designee, immediately, and investigated according to the institution's procedures. Record observations and dispositions of flasks below.

---

---

---

---

---

---

---

---

**Inspected by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

If the Site Principal Investigator, or designee, is notified of any unusual islets morphology or evidence of microbial contamination, complete the following:

**Name of Person notified:** \_\_\_\_\_

**Notified by:** \_\_\_\_\_

**Date & Time Notified:** \_\_\_\_\_

- 12.7 Post-Culture Islet Recombination – High Purity Islets
- 12.7.1 Place all the High Purity Islets T-175 culture flasks at a 45° angle and allow the islets to settle to the bottom corner for 3 to 5 minutes.
- 12.7.2 After the supernatant is observed to be clear, carefully transfer the tissue in approximately 10 mL of media from each T-175 culture flask to a T-75 flask labeled "Islets – High Purity."
- 12.7.3 Rinse the interior surfaces of each T-175 culture flask with the 20 mL of media remaining and transfer these rinses to a new T-175 flask labeled "Supernatant – High Purity."
- 12.7.4 Allow the pooled islets in the "Islets – High Purity" T-75 flask to settle for approximately 3 to 5 minutes. Remove the supernatant from the top to leave 100 mL (=100 g) of suspension in the T-75 flask. Place the supernatant into the "Supernatant – High Purity" T-175 flask.
- 12.7.5 Examine the "Supernatant – High Purity" T-175 flask under a microscope to determine if islets are present. If islets are present, transfer the supernatant to a 250 mL conical tube and centrifuge at 140 X g for 2 to 3 minutes at 2°C to 8°C. Transfer the tissue to the "Islets – High Purity" T-75 flask.

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

|   |                           |  |  |                     |
|---|---------------------------|--|--|---------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 8 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                     |

12.8 Post-Culture Islet Recombination – Middle Purity Islets

- 12.8.1 Place all the Middle Purity Islets T-175 culture flasks at a 45° angle and allow the islets to settle to the bottom corner for 3 to 5 minutes.
- 12.8.2 After the supernatant is observed to be clear, carefully transfer the tissue in approximately 10 mL of media from each T-175 culture flask to a T-75 flask labeled “Islets – Middle Purity.”
- 12.8.3 Rinse the interior surfaces of each T-175 culture flask with the 20 mL of media remaining and transfer these rinses to a new T-175 flask labeled “Supernatant – Middle Purity.”
- 12.8.4 Allow the pooled islets in the “Islets – Middle Purity” T-75 flask to settle for approximately 3 – 5 minutes. Remove the supernatant from the top to leave 100 mL (=100 g) of suspension in the T-75 flask. Place the supernatant into the “Supernatant – Middle Purity” T-175 flask.
- 12.8.5 Examine the “Supernatant – Middle Purity” T-175 flask under a microscope to determine if islets are present. If islets are present, transfer the supernatant to a 250 mL conical tube and centrifuge at 140 X g for 2 to 3 minutes at 2°C to 8°C. Transfer the tissue to the “Islets – Middle Purity” T-75 flask.

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

12.9 Post-Culture Islet Recombination – Low Purity Islets

- 12.9.1 Place all the Low Purity Islets T-175 culture flasks at a 45° angle and allow the islets to settle to the bottom corner for 3 to 5 minutes.
- 12.9.2 After the supernatant is observed to be clear, carefully transfer the tissue in approximately 10 mL of media from each T-175 culture flask to a T-75 flask labeled “Islets – Low Purity.”
- 12.9.3 Rinse the interior surfaces of each T-175 culture flask with the 20 mL of media remaining and transfer these rinses to a T-175 flask labeled “Supernatant – Low Purity.”
- 12.9.4 Allow the pooled islets in the “Islets – Low Purity” T-175 flask to settle for approximately 3 to 5 minutes. Remove the supernatant from the top to leave 100 mL (=100 g) of suspension in the T-75 flask. Place the supernatant into the “Supernatant – Low Purity” T-175 flask.
- 12.9.5 Examine the “Supernatant – Low Purity” T-175 flask under a microscope to determine if islets are present. If islets are present, transfer the supernatant to a 250 mL conical tube and centrifuge at 140 X g for 2 to 3 minutes at 2°C to 8°C. Transfer the tissue to the “Islets – Low Purity” T-75 flask.

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

|   |                           |  |  |                     |
|---|---------------------------|--|--|---------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 9 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                     |

12.10 Estimate the Settled Tissue Volume in the High, Middle and Low Purity Islets flasks

12.10.1 Allow the tissue to settle in the corner of the High Purity T-75 flask for 3 to 5 minutes.

12.10.2 Gently aspirate the tissue into a 10 mL glass pipet.

12.10.3 Allow the tissue to settle in the pipet while holding it vertically for 3 to 5 minutes.

12.10.4 Estimate the Settled Tissue Volume from the pipet and record data on the table in Section 12.12.

12.10.5 Re-suspend the tissue in the T-75 flask.

12.10.6 Repeat steps 12.10.1 to 12.10.5 for the Middle and Low Purity islets flasks.

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

12.11 Wash Tissue in Preparation for Loading into Transplant Bags

12.11.1 Allow the tissue in each T-75 flask (High, Middle and Low Purity) to settle for 3 to 5 minutes.

12.11.2 Transfer each supernatant to 250 mL conical tubes and centrifuge at 140 X g for 3 to 5 minutes.

12.11.3 Wash the settled tissue in each T-75 with approximately 100 mL CIT Transplant Wash Media.

12.11.4 Remove the supernatant from each 250 mL conical tube and return any tissue to the appropriate T-75 flask.

12.11.5 Bring the volume in each T-75 flask (High, Middle, and Low Purity) to 100 to 200 mL in CIT Transplant Media after the second wash.

12.11.6 Take a sample of each supernatant for a Gram Stain according to the institution's procedure and send it to the appropriate lab. Report the results in Section 12.12.

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

|  |                    |                                     |                                 |               |
|--|--------------------|-------------------------------------|---------------------------------|---------------|
| Document No.<br>SOP 3101, B02-2A   | Revision No.<br>04 | Effective Date<br>04 September 2009 | Supersedes Date<br>21 July 2009 | Page 10 of 32 |
| Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE) |                    |                                     |                                 |               |

12.12 The Final Product composition is based on the Settled Tissue Volume and the Gram Stain results recorded in the table, below. Determine and record which flasks to combine, if any, so that:

- If there is  $\leq 7.5$  mL Total Settled Tissue Volume, all tissue may be combined into one Final Product T-75 flask.
- There is  $\leq 7.5$  mL of Settled Tissue Volume in **any one** Final Product T-75 flask.
- There is  $\leq 15$  mL of total Settled Tissue Volume in **all** Final Product T-75 flasks.

| Purity Level | Settled Tissue Volume (mL) | Gram Stain Result* | Disposition<br>Identify which flasks will be combined or not combined for transplant, and which will be used for research or discarded. |
|--------------|----------------------------|--------------------|---|
| High         |                            |                    |   |
| Middle       |                            |                    |   |
| Low          |                            |                    |   |
| Total        |                            |                    |   |

\*The Gram Stain results are reported on the Certificates of Analysis.

Determined by: \_\_\_\_\_ Date: \_\_\_\_\_

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

If a positive Gram Stain result is reported for any purity level, immediately notify the Site Principal Investigator, or designee.

If the Site Principal Investigator, or designee, is notified of a positive Gram Stain result, complete the following:

Name of Person notified: \_\_\_\_\_

Notified by: \_\_\_\_\_

Date & Time Notified: \_\_\_\_\_

Deviation Number: \_\_\_\_\_

|  |                    |                                     |                                 |               |
|--|--------------------|-------------------------------------|---------------------------------|---------------|
| Document No.<br>SOP 3101, B02-2A   | Revision No.<br>04 | Effective Date<br>04 September 2009 | Supersedes Date<br>21 July 2009 | Page 11 of 32 |
| Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE) |                    |                                     |                                 |               |

- 12.13 Take two 100  $\mu$ L samples of each purity level and perform counts and calculations. Attach spreadsheet if used.

#### Post-culture Islets Counts

|                       | High Purity |      |     | Middle Purity |      |     | Low Purity |      |     |
|-----------------------|-------------|------|-----|---------------|------|-----|------------|------|-----|
| Sample Volume         | $\mu$ L     |      |     | $\mu$ L       |      |     | $\mu$ L    |      |     |
| Total Volume          | mL          |      |     | mL            |      |     | mL         |      |     |
| Dilution Factor       |             |      |     |               |      |     |            |      |     |
| Diameter, Factor      | Counts      | Avg. | IEQ | Counts        | Avg. | IEQ | Counts     | Avg. | IEQ |
| 50 – 100,<br>0.167    |             |      |     |               |      |     |            |      |     |
| 101 – 150,<br>0.648   |             |      |     |               |      |     |            |      |     |
| 151 – 200,<br>1.685   |             |      |     |               |      |     |            |      |     |
| 201 – 250,<br>3.500   |             |      |     |               |      |     |            |      |     |
| 251 – 300,<br>6.315   |             |      |     |               |      |     |            |      |     |
| 301 – 350,<br>10.352  |             |      |     |               |      |     |            |      |     |
| > 350, 15.833         |             |      |     |               |      |     |            |      |     |
| Total                 |             |      |     |               |      |     |            |      |     |
| % Trapped             |             |      |     |               |      |     |            |      |     |
| % Fragmented          |             |      |     |               |      |     |            |      |     |
| Purity (%)            |             |      |     |               |      |     |            |      |     |
| Islet Quality Grade*  |             |      |     |               |      |     |            |      |     |
| Technicians' Initials |             |      |     |               |      |     |            |      |     |



|  |                    |                                     |                                 |               |
|--|--------------------|-------------------------------------|---------------------------------|---------------|
| Document No.<br>SOP 3101, B02-2A   | Revision No.<br>04 | Effective Date<br>04 September 2009 | Supersedes Date<br>21 July 2009 | Page 12 of 32 |
| Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE) |                    |                                     |                                 |               |

#### Post-culture Islets Calculations

|  | High Purity | Middle Purity | Low Purity | Total |
|--|-------------|---------------|------------|-------|
| Post-culture IPN                                 |             |               |            |       |
| Post-culture IEQ                                 |             |               |            |       |
| Pre-purification IEQ<br>(Section 7.5.2)          |             |               |            |       |
| IEQ Recovery (%)<br>(from Pre-purification IEQ)  |             |               |            |       |
| Post-purification IEQ<br>(Section 10.2)          |             |               |            |       |
| IEQ Recovery (%)<br>(from Post-purification IEQ) |             |               |            |       |
| IEQ/g of trimmed pancreas<br>(Section 5.8)       |             |               |            |       |
| Comments   |             |               |            |       |

\*See Islet Quality Grade Note at the end of PBR Part 1, Section 10.2, for guidelines

Calculated by: \_\_\_\_\_ Date: \_\_\_\_\_

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

Total Post-purification Islets Count: \_\_\_\_\_ IEQ

Total Post-culture Islets Count: \_\_\_\_\_ IEQ

Percent Change: \_\_\_\_\_%

Calculated by: \_\_\_\_\_ Date: \_\_\_\_\_

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

If the Post-culture Islets Count is > 30% less than the Post-purification Islets Count, Section 10.2, notify the Site Principal Investigator, or designee, immediately.

If the Site Principal Investigator, or designee, is notified of > 30% decrease in IEQ, complete the following:

Name of Person notified: \_\_\_\_\_

Notified by: \_\_\_\_\_

Date & Time Notified: \_\_\_\_\_

|   |                           |  |  |                      |
|---|---------------------------|--|--|----------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 13 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                      |

#### 12.14 Post-culture Sampling of High Purity Islets Suspension

Based on the Post-culture count, Section 12.13, take samples of the High Purity Islets suspension according to the table below and record test results in Section 17.2, the Certificate of Analysis and Section 20.0, as required.

| SAMPLE QUANTITY         | REQUIRED FOR CERTIFICATE OF ANALYSIS                                       | SAMPLE REMOVED (mL) |
|-------------------------|--|---------------------|
| Suspension, 400 IEQ     | Post-culture Glucose Stimulated Insulin Release                            |                     |
|                         | <b>REQUIRED PRODUCT CHARACTERIZATION,<br/>FOR INFORMATION ONLY</b>         |                     |
| Suspension, 4,000 IEQ   | <i>In vivo</i> (Nude Mouse) Islets Function                                |                     |
|                         | <b>OPTIONAL PRODUCT CHARACTERIZATION,<br/>FOR INFORMATION ONLY</b>         |                     |
| Suspension, 3 X 100 IEQ | Post-culture DNA Content*  |                     |
| Suspension, 500 IEQ     | ATP/DNA  |                     |
| Suspension, 3 X 100 IEQ | Nuclei Measurement*  |                     |
| Suspension, 5,000 IEQ   | OCR/DNA*   |                     |
| Suspension, 5,000 IEQ   | Molecular Profiling*   |                     |
| Suspension, 500 IEQ     | Islets Fraction*   |                     |
|                         | Total Volume Removed<br>from High Purity Islets Suspension                 |                     |
|                         | Volume of High Purity Islets Suspension<br>Before Sampling (Section 12.13) |                     |
|                         | Remaining Volume<br>of High Purity Islets Suspension                       |                     |

\*Note: Follow instructions in the Laboratory Manual for preparation and shipment of samples.

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

12.15 Label with at least the following information one Purified Human Pancreatic Islets product infusion bag for each T-75 flask remaining after combining in Section 12.12, that will be transplanted:

- "Human Islets" or "Human Islet Product"
- Lot Number
- Donor Identification (UNOS or DDD) Number
- Donor Blood Type
- Total IEQ in Bag
- "Bag X of Y"
- Recipient Name
- Recipient Medical Record Number
- Recipient Study ID #
- Recipient Blood Type
- "Sterility testing has not been completed."
- "Biohazard: Human Tissue"

|   |                           |  |  |                      |
|---|---------------------------|--|--|----------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 14 of 32</b> |
| <b>Document Title:</b> PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE) |                           |  |  |                      |

- “New drug. Limited by law to investigational use only”
- Suspension Volume
- Name of the Manufacturing Institution
- FDA Registration Number, if available
- “BB-IND 9336”
- Storage Temperature
- “Contains Heparin, Total Units: \_\_\_\_\_”

Additional information may be added as required by the institution’s procedures.

Make three identical labels for each bag. Place one on the bag, one in the space below, or on the back of this page, if necessary, and send one with the product bag with an instruction to affix it to the recipient’s medical record chart.

**Labeled by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Checked by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

|   |                           |  |  |                      |
|---|---------------------------|--|--|----------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 15 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                      |

## 12.16 Combine the Islets Suspensions

- 12.16.1 If, according to the plan in Section 12.12, there will be one infusion bag, combine all islets into one T-75 flask rinsing the emptied flasks with CIT Transplant Media. The volume in single T-75 flask after combination should be 100 mL. Combine by settling and removing supernatant as in Section 12.11, as necessary.

Final Volume in one T-75 flask: \_\_\_\_\_ mL

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

- 12.16.2 If, according to the plan in Section 12.12, there will be two infusion bags, combine the islets into two T-75 flasks according to the plan, rinsing the emptied flasks with CIT Transplant Media. The volume in the T-75 flasks after combination should be 100 mL each. Combine by settling and removing supernatant as in Section 12.11, as necessary.

Final Volume in T-75 flask #1: \_\_\_\_\_ mL

Final Volume in T-75 flask #2: \_\_\_\_\_ mL

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

- 12.16.3 If, according to the plan in Section 12.12, there will be three infusion bags, combine the islets into three T-75 flasks according to the plan. The volume in the T-75 flasks after combination should be 100 mL each. Combine by settling and removing supernatant as in Section 12.11, as necessary.

Final Volume in T-75 flask #1: \_\_\_\_\_ mL

Final Volume in T-75 flask #2: \_\_\_\_\_ mL

Final Volume in T-75 flask #3: \_\_\_\_\_ mL

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

- 12.17 Label sample containers for the release and characterization testing samples according to the institution's procedures.

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

## 12.18 Sampling and Testing of Final Product containers

- 12.18.1 Estimate the Tissue Volume in the final product containers

- Allow the tissue to settle in the corner of T-75 flask #1 for 3 to 5 minutes.
- Gently aspirate the tissue into a sterile 10 mL glass pipet.
- Allow the tissue to settle in the pipet while holding it vertically for 3 to 5 minutes.
- Estimate the settled tissue volume from the pipet and record result in the table below.
- Re-suspend the tissue in the T-75 flask.
- Repeat these steps for other T-75 flasks.

|  |                    |                                     |                                 |               |
|--|--------------------|-------------------------------------|---------------------------------|---------------|
| Document No.<br>SOP 3101, B02-2A   | Revision No.<br>04 | Effective Date<br>04 September 2009 | Supersedes Date<br>21 July 2009 | Page 16 of 32 |
| Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE) |                    |                                     |                                 |               |

|                               | T-75 FLASK #1 | T-75 FLASK #2 | T-75 FLASK #3 |
|-------------------------------|---------------|---------------|---------------|
| SETTLED TISSUE<br>VOLUME (mL) |               |               |               |

Report these results on the Interim and Final Certificates of Analysis.

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

12.18.2 Sample the suspension(s) in the final product T-75 flask(s) [Sample the supernatant(s) for the Endotoxin test only] before filling the infusion bags, and send the samples to the appropriate laboratory for the tests indicated in the table below. Report the test results in Sections 14 and 20.0, and on the Certificates of Analysis, as indicated.

**Note:** Samples for Islets Identity and Quantity are not taken here for purity levels (High, Middle, and/or Low) that have not been combined with other purity levels for transplant. Results of the Post-culture Identity and Counts samples taken in Section 12.13, are used for the Certificates of Analysis.

| SAMPLE TYPE & QUANTITY   | TESTS                          | SAMPLE REMOVED (mL) |            |            | TESTING LAB                             |
|--|--------------------------------|---------------------|------------|------------|---|
| Interim Certificate of Analysis &<br>Certificate of Analysis                           |                                | T-75<br>#1          | T-75<br>#2 | T-75<br>#3 |   |
| 2 X 100 µL/Each Final Product<br>T-75 Flask  | Islet Identity<br>and Quantity | (A)                 |            |            |   |
| 100 IEQ/Each Final Product T-75 Flask  | Viability                      | (B)                 |            |            |   |
| 1 mL of Supernatant/Each Final Product<br>T-75 Flask                                   | Endotoxins                     | (C)                 |            |            |   |
| Certificate of Analysis  |                                |                     |            |            |   |
| 3 mL/Each Final Product<br>T-75 Flask  | Sterility<br>21 CFR 610.12     | (D)                 |            |            |   |
| Required Product Characterization,<br>For Information Only                             |                                |                     |            |            |   |
| 1,000 IEQ/Each Final Product<br>T-75 Flask   | Cell<br>Composition            | (E)                 |            |            | University of Miami*                    |
| 500 to 1,000 IEQ/Each Final Product<br>T-75 Flask                                      | MCP-1 &<br>Tissue Factor       | (F)                 |            |            | Uppsala University<br>Hospital, Sweden* |
| Optional Product Characterization,<br>For Information Only                             |                                |                     |            |            |   |
| 2,000 IEQ/Each Final Product<br>T-75 Flask   | β-cell Viability               | (G)                 |            |            |   |
| Total Volume Removed (mL) (H) = Σ (A) through (G)                                      |                                | (H)                 |            |            |   |
| Remaining Suspension Volume (mL) (J) = 100 mL – (H)<br>Required in Section 14.3, below |                                | (J)                 |            |            |   |
| Sampled by: _____  |                                | Date: _____         |            |            | Verified by: _____                      |
|  |                                |                     |            |            | Date: _____                             |

\*Note: Follow instructions in the Laboratory Manual for preparation and shipment of samples for Cell Composition, and for MCP-1 and Tissue Factor.



|   |                           |  |  |                      |
|---|---------------------------|--|--|----------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 17 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                      |

12.19 Perform counts and calculations (Portions of product that are not combined with other portions are not counted again. Their values from Section 12.13 are used.)

#### Final Product Islets Counts

|                              | Final Product T-75 Flask #1 |      |     |  | Final Product T-75 Flask #2 |      |     |  | Final Product T-75 Flask #3 |      |     |  |
|------------------------------|-----------------------------|------|-----|--|-----------------------------|------|-----|--|-----------------------------|------|-----|--|
| <b>Sample Volume</b>         | μL                          |      |     |  | μL                          |      |     |  | μL                          |      |     |  |
| <b>Total Volume</b>          | mL                          |      |     |  | mL                          |      |     |  | mL                          |      |     |  |
| <b>Dilution Factor</b>       |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| <b>Diameter, Factor</b>      | Counts                      | Avg. | IEQ |  | Counts                      | Avg. | IEQ |  | Counts                      | Avg. | IEQ |  |
| 50 – 100,<br>0.167           |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| 101 – 150,<br>0.648          |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| 151 – 200,<br>1.685          |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| 201 – 250,<br>3.500          |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| 251 – 300,<br>6.315          |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| 301 – 350,<br>10.352         |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| > 350, 15.833                |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| <b>Total</b>                 |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| <b>% Trapped</b>             |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| <b>% Fragmented</b>          |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| <b>Purity (%)</b>            |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| <b>Islet Quality Grade*</b>  |                             |      |     |  |                             |      |     |  |                             |      |     |  |
| <b>Technicians' Initials</b> |                             |      |     |  |                             |      |     |  |                             |      |     |  |

|  |                    |                                     |                                 |               |
|--|--------------------|-------------------------------------|---------------------------------|---------------|
| Document No.<br>SOP 3101, B02-2A   | Revision No.<br>04 | Effective Date<br>04 September 2009 | Supersedes Date<br>21 July 2009 | Page 18 of 32 |
| Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE) |                    |                                     |                                 |               |

## Final Product Islets Calculations

|                   | Final Product T-75 Flask<br>#1 | Final Product T-75 Flask<br>#2 | Final Product T-75 Flask<br>#3 |
|-------------------|--------------------------------|--------------------------------|--------------------------------|
| Final Product IPN |                                |                                |                                |
| Final Product IEQ |                                |                                |                                |
| Comments          |                                |                                |                                |

\*See Islets Quality Grade Note at the end of PBR Part 1, Section 10.2 for guidelines

Total Final Product IEQ: \_\_\_\_\_

Total IEQ/g of trimmed pancreas (PBR Part 1, Section 5.8): \_\_\_\_\_

Calculated by: \_\_\_\_\_ Date: \_\_\_\_\_

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

12.20 Set up the labeled product bag(s), 150 mL rinse bag(s), 60 mL syringe(s) in the BSC as follows:

- Connect the tubing from the 150 mL rinse bag to the Ricordi Infusion bag.
- Clamp off the line connecting the bags with a hemostat at both ends.
- Place a syringe in ring stand and remove its plunger.
- Connect the syringe to the Luer lock port of the Ricordi Infusion bag.
- Repeat this setup for the 2<sup>nd</sup> and 3<sup>rd</sup> bag systems, if the final tissue volume warrants multiple bags.

Performed by: \_\_\_\_\_ Date: \_\_\_\_\_

12.21 Calculation of Heparin Quantity Addition

\*\*\*\*\*

*Heparin is not a part of the product. It is added to the product at the discretion of the recipient's physician.*

\*\*\*\*\*

To the final product add 70 Units of heparin per kg of recipient body weight.

Recipient Body Weight (Section 12.3): \_\_\_\_\_ kg

Heparin Concentration: \_\_\_\_\_ units/mL

Divide the heparin equally among the infusion bags.

\_\_\_\_\_ kg X 70 U/kg/ \_\_\_\_\_ # of bags = \_\_\_\_\_ Units of Heparin to add  
to each product bag

\_\_\_\_\_ Units of Heparin to add/ \_\_\_\_\_ U/mL = \_\_\_\_\_ mL of Heparin to add  
to each product bag

Calculated by: \_\_\_\_\_ Date: \_\_\_\_\_

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

|   |                           |  |  |                      |
|---|---------------------------|--|--|----------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 19 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                      |

## 12.22 Filling Infusion and Rinse Bags #1

12.22.1 Add 100 mL of CIT Transplant Media to Infusion Bag #1. Unclamp tubing to drain the media from the infusion bag to the rinse bag. Remove all air from rinse bag and re-clamp tubing.

12.22.2 Transfer the tissue in 100 mL of CIT Transplant Media from the flask to Infusion Bag #1 through the syringe.

12.22.3 Record the time as Infusion Bag #1 Filling Start Time: \_\_\_\_\_

12.22.4 If heparin is to be added to the product, add the amount of heparin calculated in Section 12.21, to Infusion Bag #1 at this point.

Units of Heparin added to Infusion Bag #1: \_\_\_\_\_ units

Volume of Heparin added to Infusion Bag #1: \_\_\_\_\_ mL

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

12.22.5 Add 50 mL of CIT Transplant Media to the T-75 flask, rinse the surfaces of the flask with this media, and transfer this rinse media into the infusion bag.

12.22.6 Rinse the T-75 flask again with another 50 mL of CIT Transplant Media, and transfer this rinse media into the infusion bag. After transferring the entire final product to the infusion bag remove the air using a “burping” technique and clamp the port with a hemostat so that no air enters the bag.

12.22.7 Record the time as the Infusion Bag #1 Filling End Time: \_\_\_\_\_

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

## 12.23 Filling Infusion and Rinse Bags #2

12.23.1 Add 100 mL of CIT Transplant Media to Infusion Bag #2. Unclamp tubing to drain the media from the infusion bag to the rinse bag. Remove all air from rinse bag and re-clamp tubing.

12.23.2 Transfer the tissue in 100 mL of CIT Transplant Media from the flask to the Infusion Bag #2 through the syringe.

12.23.3 Record the time as Infusion Bag #2 Filling Start Time: \_\_\_\_\_

12.23.4 If heparin is to be added to the product, add the amount of heparin calculated in Section 12.21, to Infusion Bag #2 at this point.

Units of Heparin added to Infusion Bag #2: \_\_\_\_\_ units

Volume of Heparin added to Infusion Bag #2: \_\_\_\_\_ mL

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

|   |                           |  |  |                      |
|---|---------------------------|--|--|----------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 20 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                      |

12.23.5 Add 50 mL of CIT Transplant Media to the T-75 flask, rinse the surfaces of the flask with this media, and transfer this rinse media into the infusion bag.

12.23.6 Rinse the T-75 flask again with another 50 mL of CIT Transplant Media, and transfer this rinse media into the infusion bag. After transferring the entire final product to the infusion bag remove the air using a “burping” technique and clamp the port with a hemostat so that no air enters the bag.

12.23.7 Record the time as the Infusion Bag #2 Filling End Time: \_\_\_\_\_

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

#### 12.24 Filling Infusion and Rinse Bags #3

12.24.1 Add 100 mL of CIT Transplant Media to Infusion Bag #3. Unclamp tubing to drain the media from the infusion bag to the rinse bag. Remove all air from rinse bag and re-clamp tubing.

12.24.2 Transfer the tissue in 100 mL of CIT Transplant Media from the flask to Infusion Bag #3 through the syringe.

12.24.3 Record the time as Infusion Bag #3 Filling Start Time: \_\_\_\_\_

12.24.4 If heparin is to be added to the product, add the amount of heparin calculated in Section 12.21, to Infusion Bag #3 at this point.

Units of Heparin added to Infusion Bag #3: \_\_\_\_\_ units

Volume of Heparin added to Final Product Bag #3: \_\_\_\_\_ mL

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

12.24.5 Add 50 mL of CIT Transplant Media to the T-75 flask, rinse the surfaces of the flask with this media, and transfer this rinse media into the infusion bag.

12.24.6 Rinse the T-75 flask again with another 50 mL of CIT Transplant Media, and transfer this rinse media into the infusion bag. After transferring the entire final product to the infusion bag remove the air using a “burping” technique and clamp the port with a hemostat so that no air enters the bag.

12.24.7 Record the time as Infusion Bag #3 Filling End Time: \_\_\_\_\_

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

|   |                           |  |  |                      |
|---|---------------------------|--|--|----------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 21 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                      |

- 12.25 Inspect each infusion bag to ensure that it is intact, there are no leaks, the label is legible, and the contents are a light yellow to amber liquid with visible islets in each bag. These observations are reported on the Interim Certificate of Analysis and the Certificate of Analysis.

Does each product infusion bag meet these criteria?

Bag #1:            Yes                      No                      (Circle One)

Bag #2:            Yes                      No                      (Circle One)

Bag #3:            Yes                      No                      (Circle One)

If any infusion bag does not meet these criteria, the Laboratory Director, or designee, must be notified immediately, and they must initiate an investigation according to the institution's procedures. The process for reporting a deviation to the CMCMC as defined in DAIT SOP 3200 must also be followed.

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

If the Laboratory Director, or designee, is notified, complete the following:

**Name of Person notified:** \_\_\_\_\_

**Notified by:** \_\_\_\_\_

**Date & Time Notified:** \_\_\_\_\_, \_\_\_\_\_

- 12.26 Place the product infusion bags in a cooler with following:

- Absorbent material
- Room temperature pack
- Temperature monitor
- Infusion Set

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_



|  |                    |                                     |                                 |               |
|--|--------------------|-------------------------------------|---------------------------------|---------------|
| Document No.<br>SOP 3101, B02-2A   | Revision No.<br>04 | Effective Date<br>04 September 2009 | Supersedes Date<br>21 July 2009 | Page 22 of 32 |
| Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE) |                    |                                     |                                 |               |

### 13.0 CHECKLIST OF RECORDS FILED WITH THIS PRODUCTION BATCH RECORD

#### 13.1 Required Solution and Media Preparation Records

| MPBR<br>SECTION | DAIT<br>SOP 3106, | MEDIA  | PRESENT? |    |
|-----------------|-------------------|--|----------|----|
|                 |                   |  | YES      | NO |
| 5.4             | B01               | CIT Digestion Solution   |          |    |
| 5.9             | B11               | CIT Enzyme Solution  |          |    |
| 7.4.1           | B02               | CIT Purification Solution  |          |    |
| 7.4.1           | B12               | CIT Wash Solution  |          |    |
| 8.1             | B10               | CIT Purification Density Gradients   |          |    |
| 9.1             | B10               | CIT Purification Density Gradients<br>(Supplementary Purification, if performed) |          |    |
| 10.1            | B04               | CIT Culture Media  |          |    |
| 12.4.2          | B05               | CIT Transplant Wash Media  |          |    |
| 12.4.2          | B06               | CIT Transplant Media   |          |    |

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

#### 13.2 Required Lists

| MPBR<br>SECTION | LISTS   | PRESENT? |    |
|-----------------|---|----------|----|
|                 |   | YES      | NO |
| 3.1.2           | Personnel participating in this manufacturing process |          |    |
| 3.1.4           | Sterilized Items                                      |          |    |
| 3.1.5           | Equipment   |          |    |
| 3.1.6           | Disposable Items                                      |          |    |

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

#### 13.3 Required Test Reports (Results not recorded in previous Sections of this Batch Record)

| MPBR<br>SECTION | TEST REPORTS  | PRESENT? |    |
|-----------------|---|----------|----|
|                 |   | YES      | NO |
| 12.11.6         | Gram Stain  |          |    |
| 12.18.2         | Final Product Viability                               |          |    |
| 12.18.2         | Final Product Endotoxin                               |          |    |
| 12.18.2         | Pre-culture Sample Glucose Stimulated Insulin Release |          |    |

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

|  |                    |                                     |                                 |               |
|--|--------------------|-------------------------------------|---------------------------------|---------------|
| Document No.<br>SOP 3101, B02-2A   | Revision No.<br>04 | Effective Date<br>04 September 2009 | Supersedes Date<br>21 July 2009 | Page 23 of 32 |
| Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE) |                    |                                     |                                 |               |

#### 13.4 Deviation and Discrepancy Investigation Reports

Ensure that all Deviation and Discrepancy Reports related to this Batch Record are attached and approved according to the institution's procedures.

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

### 14.0 Pre-transplant Test Results

14.1 From the tests conducted on the samples from Section 12.18, enter the results in the table below.

| FINAL PRODUCT T-75 FLASKS   | #1 | #2 | #3 | TOTAL |
|---|----|----|----|-------|
| Settled Tissue Volume (mL)*   |    |    |    |       |
| Suspension Volume (mL)*   |    |    |    |       |
| Islets Identity (Yes/No)*   |    |    |    |       |
| Islets Equivalents (IEQ)  |    |    |    |       |
| Islets Quantity (IEQ/kg)*<br>(Calculate in Section 14.2, below)               |    |    |    |       |
| Islets Concentration (IEQ/mL Tissue)*<br>(Calculate in Section 14.3, below)   |    |    |    |       |
| Viability (%)*  |    |    |    |       |
| Endotoxins Concentration (EU/mL)  |    |    |    |       |
| Endotoxins<br>(EU/kg Recipient Weight)*<br>(Calculate in Section 14.4, below) |    |    |    |       |

\*These results are also reported on the Interim and Final Certificates of Analysis.

14.2 From the Islets Equivalents in Section 14.1, above, and the Recipient Body Weight (kg) in Section 12.3, above, calculate the Islets Quantity (IEQ/kg) in each T-75 Flask and their sum, and record the results in the table above:

$\frac{\text{Islets Equivalents (IEQ)}}{\text{Recipient Body Weight (kg)}} = \text{Islets Quantity (IEQ/kg)}$

| FINAL PRODUCT T-75 FLASKS | Islets Equivalents (IEQ) | Recipient body Weight (kg) | Islets Quantity (IEQ/kg) |
|---------------------------|--------------------------|----------------------------|--------------------------|
| 1                         |                          |                            |                          |
| 2                         |                          |                            |                          |
| 3                         |                          |                            |                          |

Entered and calculated by: \_\_\_\_\_

Date: \_\_\_\_\_

Verified by: \_\_\_\_\_

Date: \_\_\_\_\_

|   |                           |  |  |                      |
|---|---------------------------|--|--|----------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 24 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                      |

- 14.3 From the Islets Equivalents and the Settled Tissue Volumes in Section 14.1, above, calculate the Islets Concentration in each T-75 Flask and their sum, and record the results in the table above:

$$\frac{\text{Islets Equivalents (IEQ)}}{\text{Settled Tissue Volume (mL)}} = \text{Islets Concentration (IEQ/mL Tissue)}$$

| FINAL PRODUCT<br>T-75 FLASKS | Islets Equivalents<br>(IEQ) | Settled Tissue Volume<br>(mL) | Islets Concentration<br>(IEQ/mL) |
|------------------------------|-----------------------------|-------------------------------|----------------------------------|
| 1                            |                             |                               |                                  |
| 2                            |                             |                               |                                  |
| 3                            |                             |                               |                                  |

Entered and calculated by: \_\_\_\_\_ Date: \_\_\_\_\_

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

- 14.4 From the Endotoxins Concentration (EU/mL) in Section 14.1, the Suspension Volume (mL) in Section 12.18.2 (J), and the Recipient Body Weight (kg) in Section 12.3, above, calculate the Endotoxins (Units per kg of recipient body weight) in each T-75 Flask and the Total Endotoxins (Units per kg of recipient body weight), and record the results in the table above:

$$\frac{\text{Endotoxins Concentration (EU/mL)} \times \text{Suspension Volume (mL)}}{\text{Recipient Body Weight (kg)}} = \text{EU/kg Recipient Weight}$$

| FINAL PRODUCT<br>T-75 FLASKS | Endotoxins<br>Concentration (EU/mL) | Suspension<br>Volume (mL) | Recipient Body<br>Weight (kg) | EU/kg |
|------------------------------|-------------------------------------|---------------------------|-------------------------------|-------|
| 1                            |                                     |                           |                               |       |
| 2                            |                                     |                           |                               |       |
| 3                            |                                     |                           |                               |       |

Entered and calculated by: \_\_\_\_\_ Date: \_\_\_\_\_

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

- 14.5 Glucose Stimulated Insulin Release Test Results

| HIGH PURITY LEVEL                    | INSULIN CONCENTRATIONS |              |                   |
|--------------------------------------|------------------------|--------------|-------------------|
|                                      | LOW GLUCOSE            | HIGH GLUCOSE | STIMULATION INDEX |
| PRE-CULTURE SAMPLE<br>(SECTION 11.1) |                        |              |                   |

Report this result on the Interim Certificate of Analysis.

Recorded by: \_\_\_\_\_ Date: \_\_\_\_\_

Verified by: \_\_\_\_\_ Date: \_\_\_\_\_

|   |                           |  |  |                      |
|---|---------------------------|--|--|----------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 25 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                      |

## 15.0 PRE-TRANSPLANT BATCH RECORD REVIEW AND INTERIM APPROVAL

After the completion of all activities and records of this manufacturing process to this point, and before transplant of this batch of islets, a qualified technician, and the Laboratory Director, Operations Manager, or designee, must review the Production Batch Record (both Part 1 and Part 2A) to verify that it is complete and accurate to this point.

We have reviewed the Production Batch Record (both Part 1 and Part 2A) and verified that it is complete and accurate to this point.

\_\_\_\_\_  
Qualified Technician

Date: \_\_\_\_\_

\_\_\_\_\_  
Laboratory Director, Operations Manager, or designee

Date: \_\_\_\_\_

|   |                           |  |  |                      |
|---|---------------------------|--|--|----------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 26 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                      |

## 16.0 ISLET PRODUCT CUSTODY TRANSFER

16.1 Notify the clinical team that the islets are ready for transplant.

**Notified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_ **Time:** \_\_\_\_\_

16.2 Custody Transfer Record

File the original or a copy of the institution's product custody transfer record with this production batch record.

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

16.3 Review the product bag label(s) with a clinical team member to assure that the intended recipient and the UNOS or DDD Number are correctly identified (See Section 12.3). Report this identity verification on the Interim and Final Certificates of Analysis.

UNOS or DDD Number Correct? Yes No (Circle One)

Recipient Identity Correct? Yes No (Circle One)

**Performed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

## 17.0 POST-TRANSPLANT TEST RESULTS & REPORTS

17.1 Sterility Test Results

17.1.1 Record the 24-hour and final test results of the 21 CFR 610.12 sterility test and fungal culture on the Preservation Solution (Section 5.1) in the table below, when available.

| PRESERVATION SOLUTION | 24-HOUR RESULT | FINAL RESULT |
|-----------------------|----------------|--------------|
|                       |                |              |

If there is a positive result, record the identity of the organism(s): \_\_\_\_\_

\_\_\_\_\_  
\_\_\_\_\_

**Recorded by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_



|   |                           |  |  |                      |
|---|---------------------------|--|--|----------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 27 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                      |

- 17.1.2 Record the Final Results of the sterility test (21 CFR 610.12) and fungal culture on the samples from the Final Product T-75 Flasks (taken at Section 12.18.2) in the table below. Report these results on the Final Certificate of Analysis, when available.

| <b>FINAL PRODUCT<br/>T-75 FLASKS</b> | <b>24-HOUR RESULT</b> | <b>FINAL RESULT</b> |
|--------------------------------------|-----------------------|---------------------|
| #1                                   |                       |                     |
| #2                                   |                       |                     |
| #3                                   |                       |                     |

If there is a positive result, record the identity of the organism(s): \_\_\_\_\_

\_\_\_\_\_  
\_\_\_\_\_

**Recorded by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

If any positive result is reported, immediately notify the attending physician.

Name of Physician Notified: \_\_\_\_\_

**Notified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_ **Time:** \_\_\_\_\_

#### 17.2 Glucose Stimulated Insulin Release Test Results

| <b>HIGH PURITY LEVEL</b>                       | <b>INSULIN CONCENTRATIONS</b> |                     |                          |
|--|-------------------------------|---------------------|--------------------------|
|  | <b>LOW GLUCOSE</b>            | <b>HIGH GLUCOSE</b> | <b>STIMULATION INDEX</b> |
| <b>POST-CULTURE SAMPLE<br/>(SECTION 12.14)</b> |                               |                     |                          |

Report this result on the Certificate of Analysis.

**Recorded by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

#### 17.3 Required Test Reports (Results not recorded in previous Sections of this Batch Record)

| <b>MPBR<br/>SECTION</b> | <b>TEST REPORTS</b>                              | <b>PRESENT?</b> |           |
|-------------------------|--|-----------------|-----------|
|                         |  | <b>YES</b>      | <b>NO</b> |
| 5.1                     | Preservation Solution Sterility                  |                 |           |
| 12.14                   | Final Product Glucose Stimulated Insulin Release |                 |           |
| 12.18.2                 | Final Product Sterility                          |                 |           |

**Verified by:** \_\_\_\_\_ **Date:** \_\_\_\_\_

## 18.0 PRODUCT DISPOSITION

Was this product transplanted?                      Yes                      No                      (Circle one)

If this product was transplanted, give the Patient Study ID #: \_\_\_\_\_

If this product, or any portion of it, was not transplanted, explain why not and state its final disposition.

[illegible]

Recorded by: \_\_\_\_\_ Date: \_\_\_\_\_

## 19.0 POST-TRANSPLANT BATCH RECORD REVIEW AND FINAL APPROVAL

After completion of Sections 16, 17, and 18, above, a qualified technician, and the Laboratory Director, Operations Manager, or designee review these Sections to assure that they are complete and accurate.

We have reviewed Sections 16, 17, and 18, above, and verified that they are complete and accurate.

\_\_\_\_\_  
Qualified Technician

Date: \_\_\_\_\_

\_\_\_\_\_  
Laboratory Director, Operations Manager or designee

A qualified representative of the institution's Quality Unit must review the entire Production Batch Record (both Part 1 and Part 2A) and verify that it is complete and accurate

I have reviewed this entire Production Batch Record (both Part 1 and Part 2A) and verified that it is complete and accurate.

\_\_\_\_\_  
Quality Unit Representative

Date: \_\_\_\_\_

|  |                    |                                     |                                 |               |
|--|--------------------|-------------------------------------|---------------------------------|---------------|
| Document No.<br>SOP 3101, B02-2A   | Revision No.<br>04 | Effective Date<br>04 September 2009 | Supersedes Date<br>21 July 2009 | Page 29 of 32 |
| Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE) |                    |                                     |                                 |               |

20.0 Product Characterization Test Results (For Information Only)

Record results of the following tests in the table below. File copies of the raw data with this PBR. "FPTF" means Final Product T-75 Flask.

| SAMPLES FROM MPBR SECTION | REQUIRED TEST   | RESULT  |
|---------------------------|---|---|
| 5.8                       | Pancreas Biopsy<br>MCP-1  |   |
| 5.8                       | Pancreas Biopsy<br>Tissue Factor  |   |
| 12.14                     | <i>In Vivo</i> Islet Function<br>(Nude Mouse Assay)                       | High Purity Islets: _____<br>(Hyperglycemia Reversed, or Not Reversed)  |
| 12.18.2                   | Cell Composition<br>(Laser Scanning<br>Cytometry &<br>Immunofluorescence) | FPTF #1, $\beta$ -cells: _____ %<br>$\gamma$ -cells: _____ %<br>$\alpha$ -cells: _____ %<br>PP-cells: _____ %<br>FPTF #2, $\beta$ -cells: _____ %<br>$\gamma$ -cells: _____ %<br>$\alpha$ -cells: _____ %<br>PP-cells: _____ %<br>FPTF #3, $\beta$ -cells: _____ %<br>$\gamma$ -cells: _____ %<br>$\alpha$ -cells: _____ %<br>PP-cells: _____ % |
| 12.18.2                   | Final Product<br>MCP-1  | FPTF 1: _____<br>FPTF 2: _____<br>FPTF 3: _____   |
| 12.18.2                   | Final Product<br>Tissue Factor  | FPTF 1: _____<br>FPTF 2: _____<br>FPTF 3: _____   |
| SAMPLES FROM MPBR SECTION | OPTIONAL TEST   | RESULT  |
| 11.1                      | Pre-culture<br>DNA Content  | High Purity Islets: _____ $\mu$ g DNA   |
| 11.1                      | Pre-culture<br>Nuclei Measurement   | _____ nuclei  |
| 12.14                     | Post-culture<br>DNA Content   | High Purity Islets: _____ $\mu$ g DNA   |
| 12.14                     | Post-culture<br>Nuclei Measurement  | _____ nuclei  |
| 12.14                     | ATP/DNA Ratio   |   |
| 12.14                     | OCR/DNA   | _____ nmol O <sub>2</sub> /min/mg DNA   |
| 12.14                     | Molecular Profiling   |   |
| 12.14                     | Islet Fraction  |   |
| 12.18.2                   | $\beta$ -Cell Viability<br>(Flow Cytometry)                               | FPTF #1: _____ %<br>FPTF #2: _____ %<br>FPTF #3: _____ %  |

Recorded by: \_\_\_\_\_

Date: \_\_\_\_\_


Verified by: \_\_\_\_\_

Date: \_\_\_\_\_

|  |                    |                                     |                                 |               |
|--|--------------------|-------------------------------------|---------------------------------|---------------|
| Document No.<br>SOP 3101, B02-2A   | Revision No.<br>04 | Effective Date<br>04 September 2009 | Supersedes Date<br>21 July 2009 | Page 30 of 32 |
| Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE) |                    |                                     |                                 |               |

## HEMODILUTION FLOWCHART

### DONOR SPECIMEN SUITABILITY FOR INFECTIOUS DISEASE TESTING FLOWCHART



#### Definitions:

1. Blood or blood component: any part of a single-donor unit of blood separated by physical or mechanical means.
2. Colloid: a protein or polysaccharide solution that can be used to increase or maintain osmotic (oncotic) pressure in the intravascular compartment such as albumin, dextran, hetastarch; or certain blood components, such as plasma or platelets.
3. Crystalloid: a balanced salt and/or glucose solution used for electrolyte replacement or to increase intravascular volume such as saline, Ringer's lactate solution, or 5% dextrose in water.

|   |                           |  |  |                      |
|---|---------------------------|--|--|----------------------|
| <b>Document No.</b><br>SOP 3101, B02-2A   | <b>Revision No.</b><br>04 | <b>Effective Date</b><br>04 September 2009 | <b>Supersedes Date</b><br>21 July 2009 | <b>Page 31 of 32</b> |
| <b>Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE)</b> |                           |  |  |                      |

## HEMODILUTION WORKSHEET

**Instructions:** Use this worksheet when (1) no pre-transfusion sample is available and (2) the determination needs to be made if the post-transfusion sample is suitable for infectious disease testing due to transfusion or infusion.

Donor UNOS # \_\_\_\_\_ Date: \_\_\_\_\_

|   |   |
|---|---|
| Date and Time of Sampling   | a.m.      p.m.  |
| Donor Weight (kg)   | kg  |
| Plasma Volume (PV)  | Donor weight (kg): _____/0.025 = _____ mL   |
| Blood Volume (BV)   | Donor weight (kg): _____/ 0.015 = _____ mL  |
| <b>A. Total Volume of Blood transfused/48 hours</b><br><br>1 unit packed red cells = 250 mL<br><br>Date and Time of Transfusion   | RBC's transfused/48 hrs: _____ mL<br><br>Whole blood transfused / 48 hrs: _____ mL<br><br>Reconstituted blood transfusion: _____ mL<br><br><b>Total of A: _____ mL</b>  |
| <b>B. Total Volume of colloid transfused/48 hours</b><br><br>1 unit FFP = 250 mL<br>1 unit platelet pheresis = 225 mL<br>1 platelet pool = 300 mL<br><br>Date and Time of Transfusion | Dextran / 48 hrs: _____ mL<br><br>Plasma / 48 hrs: _____ mL<br><br>Platelets / 48 hrs: _____ mL<br><br>Albumin / 48 hrs: _____ mL<br><br>Hetastarch / 48 hrs: _____ mL<br><br>Other ( _____ ): _____ mL<br><br>Other ( _____ ): _____ mL<br><br><b>Total of B: _____ mL</b> |
| <b>C. Total Volume of crystalloid transfused/1 hour</b>   | Saline: _____ mL<br><br>Dextrose in Water: _____ mL<br><br>Ringer's Lactate: _____ mL<br><br>Other ( _____ ): _____ mL<br><br>Other ( _____ ): _____ mL<br><br><b>Total of C: _____ mL</b>  |



|  |                    |                                     |                                 |               |
|--|--------------------|-------------------------------------|---------------------------------|---------------|
| Document No.<br>SOP 3101, B02-2A   | Revision No.<br>04 | Effective Date<br>04 September 2009 | Supersedes Date<br>21 July 2009 | Page 32 of 32 |
| Document Title: PHPI, MPBR, PART 2A (PRODUCT CODE PHPI-A-01, ISLETS ALONE) |                    |                                     |                                 |               |

## HEMODILUTION WORKSHEET (CONTINUED)

|  |  |
|--|--|
| <p>D. Determination of Suitability</p> <p>B _____ mL + C _____ mL = _____ mL</p> <p>A _____ mL + B _____ mL + C _____ mL</p> <p>= _____ mL</p> | <p>1. Is B + C &gt; PV? (circle one)      Yes      No</p> <p>2. Is A + B + C &gt; BV? (circle one)      Yes      No</p> <p><i>If the answers to both 1 and 2 are NO, then test sample.</i></p> <p><i>If the answer to either 1 or 2 is YES, then reject donor.</i></p> |
|--|--|

Test blood sample? (circle one)              Yes              No

Donor Suitable? (circle one)              Yes              No

Recorded by : \_\_\_\_\_              Date: \_\_\_\_\_

Reviewed by : \_\_\_\_\_              Date: \_\_\_\_\_