



### Media Preparation Form

<b>Bone Marrow Processing &amp; Resuspension Media</b>				
Date / Time of Preparation:		Expiration Date:		
HBM #:		Total number of bottles prepared:		
Prepared by:				
Ingredients	Volume to be Used	Manufacturer	Expiration Date	Actual Volume
RPMI-1640 (Ca <sup>2+</sup> /Mg <sup>2+</sup> free, no Phenol Red)	1 Liter			
Human Serum Albumin (25%)	100 ml (2.5% final conc.)			
Heparin (1,000 U/ml)	10 ml			
Gentamycin (40mg/ml)	10 ml			
DNase	<input type="checkbox"/> 1 vial = 2.5ml <input type="checkbox"/> 1/3 vial = 1.66ml			
<b>IMPORTANT REMINDER:</b>				
1. Check that the pH of all media is 7.3 ± 0.1 _____ 2. Record the lot numbers of these ingredients on SOP HCPF- 001, Attachment I				
Bottles labeled by:		Filtered by (if required):		
Comments:				

**Reviewed by:** \_\_\_\_\_ **Date:** \_\_\_\_\_